

Anti-Monoacylglycerol Lipase/MGL antibody [EPR6339] - BSA and Azide free ab248009

リコンビナント RabMAb

画像数 3

製品の概要

製品名	Anti-Monoacylglycerol Lipase/MGL antibody [EPR6339] - BSA and Azide free
製品の詳細	Rabbit monoclonal [EPR6339] to Monoacylglycerol Lipase/MGL - BSA and Azide free
由来種	Rabbit
アプリケーション	適用あり: WB 適用なし: ICC/IF, IHC-P or IP
種交差性	交差種: Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
特記事項	<p>ab248009 is the carrier-free version of ab124796.</p> <p>Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>This product is compatible with the Maxpar[®] Antibody Labeling Kit from Fluidigm, without the need for antibody preparation. Maxpar[®] is a trademark of Fluidigm Canada Inc.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> <p>Mouse, Rat: We have preliminary internal testing data to indicate this antibody may not react with</p>

these species. Please contact us for more information.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	pH: 7.2 Constituent: PBS
キャリア・フリー	はい
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR6339
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab248009の使用に適用されません
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Predicted molecular weight: 33 kDa.

追加情報 Is unsuitable for ICC/IF, IHC-P or IP.

ターゲット情報

機能	Converts monoacylglycerides to free fatty acids and glycerol. Hydrolyzes the endocannabinoid 2-arachidonoylglycerol, and thereby contributes to the regulation of endocannabinoid signaling, nociception and perception of pain (By similarity). Regulates the levels of fatty acids that serve as signaling molecules and promote cancer cell migration, invasion and tumor growth.
組織特異性	Detected in adipose tissue, lung, liver, kidney, brain and heart.
パスウェイ	Glycerolipid metabolism; triacylglycerol degradation.
配列類似性	Belongs to the AB hydrolase superfamily. Monoacylglycerol lipase family.

画像



Western blot - Anti-Monoacylglycerol Lipase/MGL antibody [EPR6339] - BSA and Azide free (ab248009)


Anti-Monoacylglycerol Lipase/MGL antibody [EPR6339] ([ab124796](#)) at 1/1000 dilution + Human heart lysate at 10 µg

Secondary

HRP labelled goat anti-rabbit at 1/2000 dilution

Predicted band size: 33 kDa

This data was developed using [ab124796](#), the same antibody clone in a different buffer formulation.



Western blot - Anti-Monoacylglycerol Lipase/MGL antibody [EPR6339] - BSA and Azide free (ab248009)

Anti-Monoacylglycerol Lipase/MGL antibody [EPR6339] ([ab124796](#)) at 1/1000 dilution + Human adipose tissue lysate at 10 µg

Secondary

HRP labelled goat anti-rabbit at 1/2000 dilution

Predicted band size: 33 kDa

This data was developed using [ab124796](#), the same antibody clone in a different buffer formulation.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-Monoacylglycerol Lipase/MGL antibody
[EPR6339] - BSA and Azide free (ab248009)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors