

Mouse IgG1, kappa monoclonal [MG1-45] - Isotype control ab18448

9 References

製品の概要

製品名	Mouse IgG1, kappa monoclonal [MG1-45] - Isotype control
アプリケーション	適用あり: IHC-P, IHC-Fr, Flow Cyt, WB, IP, ELISA
免疫原	Chemical/ Small Molecule (Trinitrophenol) conjugated to KLH.
特記事項	

The ab18448, MG1-45 immunoglobulin reacts with KLH and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse and human tissues.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein G purified
特記事項(精製)	The immunoglobulin was purified by affinity chromatography.
アイソタイプ・コントロール 備考	The ab18448, MG1-45 immunoglobulin reacts with KLH and was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse and human tissues.
ポリクローナル	モノクローナル
クローン名	MG1-45
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab18448の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. Suitable for surface or intracellular staining.
WB		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration. PubMed: 18006877

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors