

Mouse IgG3, monoclonal [PPV-07] - Isotype Control ab91537

製品の概要

製品名	Mouse IgG3, monoclonal [PPV-07] - Isotype Control
特異性	ab91537 reacts with an undefined epitope on a plant pathogen. No reactivity with cell surface antigens on Human, Porcine, Mouse or Rat. It can be used as a negative control for immunofluorescence staining (e.g. Flow Cytometry) when using Mouse monoclonal antibodies to Human cell surface antigens of isotype IgG1. ab91537 enables an estimation of non specific binding of Mouse monoclonal antibodies to cell surface components in peripheral blood and tissue. It is suitable for whole blood, Ficoll separated preparations, frozen and paraffin embedded sections.
アプリケーション	適用あり: Flow Cyt
免疫原	not applicable
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: 1% BSA
精製度	Immunogen affinity purified
特記事項(精製)	Purified (pass through 0.22 µm filter).
ポリ/モノ	モノクローナル
クローン名	PPV-07
アイソタイプ	IgG3

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab91537の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		

追加情報

Flow Cyt: Use 10 µl for 10⁶ cells.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors