


Anti-Myelin Basic Protein antibody [IGX3421] ab209328

リコンビナント

★★★★★ **9 Abreviews** **8 References** 画像数 9

製品の概要

製品名	Anti-Myelin Basic Protein antibody [IGX3421]
製品の詳細	Human monoclonal [IGX3421] to Myelin Basic Protein
由来種	Human
アプリケーション	適用あり: ELISA, WB, ICC/IF, IHC-P
種交差性	交差種: Mouse, Rat, Human, Recombinant fragment 交差が予測される動物種: Horse, Cow, Cat, Dog, Pig, Chimpanzee, Macaque monkey 
免疫原	Full length native protein (purified). This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Human, mouse and rat brain tissue lysate. Myelin Basic Protein (recombinant protein) IHC-P: Mouse, rat and human brain tissue (Hippocampus) ICC/IF: SK-N-SH cells
特記事項	This product was made using <u>synthetic libraries and phage display technology</u> . This antibody is a recombinant antibody. Human monoclonal antibody.

Example of usage (reference):

Spatiotemporal Dynamics of Molecular Pathology in Amyotrophic Lateral Sclerosis

Silas Maniatis, Tarmo Aijo, Sanja Vickovic, Catherine Braine, Kristy Kang, Annelie Mollbrink, Zaneta Andrusivova, Sami Saarenpaa, Gonzalo Saiz-Castro, Miguel Cuevas, Aaron Watters, Joakim Lundberg, Richard Bonneau, Hemali Phatnani
doi: <https://doi.org/10.1101/389270>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 0.05% BSA, 40% Glycerol (glycerin, glycerine)
精製度	Protein A purified

ポリ/モノ	モノクローナル
クローン名	IGX3421
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab209328の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
WB	★★★★★ (2)	Use a concentration of 0.25 - 1 µg/ml. Detects a band of approximately 20,17 kDa (predicted molecular weight: 33 kDa).
ICC/IF		Use a concentration of 5 µg/ml. This product gave a positive signal in SKNSH cells fixed with 4% formaldehyde
IHC-P	★★★★★ (6)	Use a concentration of 0.5 - 1 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

機能	The classic group of MBP isoforms (isoform 4-isoform 14) are with PLP the most abundant protein components of the myelin membrane in the CNS. They have a role in both its formation and stabilization. The smaller isoforms might have an important role in remyelination of denuded axons in multiple sclerosis. The non-classic group of MBP isoforms (isoform 1-isoform 3/Golli-MBPs) may preferentially have a role in the early developing brain long before myelination, maybe as components of transcriptional complexes, and may also be involved in signaling pathways in T-cells and neural cells. Differential splicing events combined with optional post-translational modifications give a wide spectrum of isomers, with each of them potentially having a specialized function. Induces T-cell proliferation.
組織特異性	MBP isoforms are found in both the central and the peripheral nervous system, whereas Golli-MBP isoforms are expressed in fetal thymus, spleen and spinal cord, as well as in cell lines derived from the immune system.
関連疾患	Note=The reduction in the surface charge of citrullinated and/or methylated MBP could result in a weakened attachment to the myelin membrane. This mechanism could be operative in demyelinating diseases such as chronic multiple sclerosis (MS), and fulminating MS (Marburg disease).
配列類似性	Belongs to the myelin basic protein family.
発生段階	Expression begins abruptly in 14-16 week old fetuses. Even smaller isoforms seem to be produced during embryogenesis; some of these persisting in the adult. Isoform 4 expression is more evident at 16 weeks and its relative proportion declines thereafter.
翻訳後修飾	Several charge isomers of MBP; C1 (the most cationic, least modified, and most abundant form), C2, C3, C4, C5, C6, C7, C8-A and C8-B (the least cationic form); are produced as a result of


optional PTM, such as phosphorylation, deamidation of glutamine or asparagine, arginine citrullination and methylation. C8-A and C8-B contain each two mass isoforms termed C8-A(H), C8-A(L), C8-B(H) and C8-B(L), (H) standing for higher and (L) for lower molecular weight. C3, C4 and C5 are phosphorylated. The ratio of methylated arginine residues decreases during aging, making the protein more cationic.

The N-terminal alanine is acetylated (isoform 3, isoform 4, isoform 5 and isoform 6). Arg-241 was found to be 6% monomethylated and 60% symmetrically dimethylated.

細胞内局在

Myelin membrane. Cytoplasmic side of myelin.


画像



IHC image of Myelin Basic Protein staining in a section of formalin-fixed paraffin-embedded normal rat brain and normal rat pancreas, performed on a Leica BOND™. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6, epitope retrieval solution 1) for 20 minutes. The section was then incubated with ab209328, 1/1000 dilution, for 15 minutes at room temperature.


An HRP-conjugated goat anti-Human IgG secondary was used for 15 minutes at room temperature. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX. The inset secondary-only control image is taken from an identical assay without primary antibody.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.



ab209328 staining Myelin Basic Protein in SK-N-SH (Human neuroblastoma cell line) cells. The cells were fixed with 4% formaldehyde (10 minutes), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal donkey serum/0.3M glycine in 0.1% PBS-Tween for 1 hour. The cells were then incubated overnight at +4°C with ab209328 at a 5 µg/ml concentration, then detected with a donkey anti-human (Alexa Fluor® 488) secondary antibody at a 1/2000 dilution (shown in green). Nuclear DNA was labelled with DAPI (shown in blue), and **ab195889**, Mouse monoclonal to alpha Tubulin (Alexa Fluor® 594), at a 1/250 dilution (shown in red).


Image was taken with a confocal microscope (Leica-Microsystems, TCS SP8).



Immunocytochemistry/ Immunofluorescence - Anti-Myelin Basic Protein antibody [IGX3421] (ab209328)

ab209328 staining Myelin Basic Protein in SHSY5Y cells. The cells were fixed with 100% methanol (5 minutes), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal donkey serum/0.3M glycine in 0.1% PBS-Tween for 1 hour. The cells were then incubated overnight at +4°C with ab209328 at a 5 µg/ml concentration, then detected with a donkey anti-human (Alexa Fluor® 488) secondary antibody at a 1/1000 dilution (shown in green). Nuclear DNA was labelled with DAPI (shown in blue), and **ab195884**, Rat monoclonal to alpha Tubulin (Alexa Fluor® 647), at a 1/250 dilution (shown in red).

Image was taken with a confocal microscope (Leica-Microsystems, TCS SP8).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myelin Basic Protein antibody [IGX3421] (ab209328)

IHC image of Myelin Basic Protein staining in a section of formalin-fixed paraffin-embedded normal human hippocampus and normal human pancreas*, performed on a Leica BOND™. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6, epitope retrieval solution 1) for 20 minutes. The section was then incubated with ab209328, 1/1000 dilution, for 15 minutes at room temperature.

An HRP-conjugated goat anti-Human IgG secondary was used for 15 minutes at room temperature. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX. The inset secondary-only control image is taken from an identical assay without primary antibody.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.

*Human pancreas tissue was obtained from the Human Research Tissue Bank, supported by the NIHR Cambridge Biomedical Research Centre




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myelin Basic Protein antibody [IGX3421] (ab209328)

IHC image of Myelin Basic Protein staining in a section of formalin-fixed paraffin-embedded normal mouse brain and normal mouse pancreas, performed on a Leica BOND™. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6, epitope retrieval solution 1) for 20 minutes. The section was then incubated with ab209328, 1/1000 dilution, for 15 minutes at room temperature.

An HRP-conjugated goat anti-Human IgG secondary was used for 15 minutes at room temperature. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX. The inset secondary-only control image is taken from an identical assay without primary antibody.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.



Western blot - Anti-Myelin Basic Protein antibody [IGX3421] (ab209328)

All lanes : Anti-Myelin Basic Protein antibody [IGX3421] (ab209328) at 0.25 µg/ml

Lane 1 : Human brain tissue lysate - total protein (**ab29466**) at 10 µg

Lane 2 : Mouse brain tissue lysate at 10 µg

Lane 3 : Rat brain tissue lysate at 10 µg

Lane 4 : Myelin Basic Protein (Recombinant protein) at 0.1 µg

Secondary

All lanes : HRP conjugated Goat Anti-Human IgG (H+L) at 1/10000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 33 kDa

Observed band size: 18,23,24 kDa


Exposure time :

Lane 1 : 30 seconds.

Lanes 2-3 : 2 minutes.

Lane 4 : 8 minutes.

This blot was produced using a 4-12% Bis-tris gel under the MES buffer system. The gel was run at 200V for 35 minutes before being transferred onto a Nitrocellulose membrane at 30V for 70 minutes. The membrane was then blocked for an hour using 3% milk before being incubated with ab209328 overnight at 4°C. Antibody binding was visualised using ECL development solution **ab133406**.



ELISA - Anti-Myelin Basic Protein antibody
[IGX3421] (ab209328)

ELISA using ab209328 for 16 hours at 4°C. **ab7153** goat anti human was used as a secondary at a 1/5000 dilution for 1 hour at Room Temperature.



Western blot - Anti-Myelin Basic Protein antibody
[IGX3421] (ab209328)

Anti-Myelin Basic Protein antibody [IGX3421] (ab209328) at 1 µg/ml + Mouse brain tissue lysate at 10 µg

Secondary

HRP conjugated Goat Anti-Human IgG (H+L) at 1/10000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 33 kDa

Observed band size: 20, 17 kDa

Exposure time: 2 minutes

This blot was produced using a 4-12% Bis-tris gel under the MES buffer system. The gel was run at 200V for 35 minutes before being transferred onto a Nitrocellulose membrane at 30V for 70 minutes. The membrane was then blocked for an hour using 3% milk before being incubated with ab209328 overnight at 4°C. Antibody binding was visualised using ECL development solution **ab133406**.

Why choose a recombinant antibody?

- Research with confidence**
Consistent and reproducible results
- Long-term and scalable supply**
Recombinant technology
- Success from the first experiment**
Confirmed specificity
- Ethical standards compliant**
Animal-free production

Anti-Myelin Basic Protein antibody [IGX3421]
(ab209328)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors