

Anti-Myoglobin antibody [EP3081Y] ab77232

リコンビナント **RabMAb**


★★★★★ **1 Abreviews** **29 References** 画像数 **14**

製品の概要

製品名	Anti-Myoglobin antibody [EP3081Y]
製品の詳細	Rabbit monoclonal [EP3081Y] to Myoglobin
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P, ICC/IF
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide within Human Myoglobin aa 100-200 (C terminal). The exact sequence is proprietary. Database link: P02144
ポジティブ・コントロール	Human, rat and mouse muscle tissue lysates; human heart muscle tissue.
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb® patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
解離定数 (K _D 値)	K _D = 1.78 x 10 ⁻¹⁰ M



[Learn more about K_D](#)

バッファー	pH: 7.20
-------	----------

	Preservative: 0.01% Sodium azide
	Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.5% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EP3081Y
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab77232の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	1/1000 - 1/5000. Predicted molecular weight: 17 kDa.
IHC-P		1/250 - 1/500. Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol. Use of an HRP/AP polymerized secondary antibody is recommended. See <u>IHC antigen retrieval protocols</u> .
ICC/IF		Use at an assay dependent concentration.

ターゲット情報

機能	Serves as a reserve supply of oxygen and facilitates the movement of oxygen within muscles.
配列類似性	Belongs to the globin family.

画像



All lanes : Anti-Myoglobin antibody [EP3081Y] (ab77232) at 1/5000 dilution (Purified)

Lane 1 : Human muscle lysate with 5% NFDM/TBST

Lane 2 : Mouse muscle lysate with 5% NFDM/TBST

Lane 3 : Rat muscle lysate with 5% NFDM/TBST

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution (Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated)


Predicted band size: 17 kDa

Observed band size: 17 kDa

Exposure time: 1 second




Purified ab77232 staining Myoglobin in Human cardiac muscle tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with paraformaldehyde; antigen retrieval was by heat mediation in a Tris/EDTA buffer (**ab93684** Tris/EDTA buffer, pH 9.0). Samples were incubated with primary antibody (1/1000). A HRP-conjugated Goat Anti-Rabbit IgG H&L (ready to use) was used as the secondary antibody. Counterstained with Hematoxylin.




Immunocytochemistry/ Immunofluorescence - Anti-Myoglobin antibody [EP3081Y] (ab77232)

Purified ab77232 staining Myoglobin in the C2C12 (mouse myoblasts myoblast) cell line by ICC/IF (Immunocytochemistry/immunofluorescence). Cells were fixed with methanol. Samples were incubated with primary antibody and an Alexa Fluor® 488-conjugated Goat anti-rabbit IgG ([ab150077](#)) was used as the secondary antibody. Counterstained with [ab195889](#) Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) and DAPI.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myoglobin antibody [EP3081Y] (ab77232)

Purified ab77232 staining Myoglobin in Mouse cardiac muscle tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with paraformaldehyde; antigen retrieval was by heat mediation in a Tris/EDTA buffer ([ab93684](#) Tris/EDTA buffer, pH 9.0). Samples were incubated with primary antibody (1/1000). A HRP-conjugated Goat Anti-Rabbit IgG H&L (ready to use) was used as the secondary antibody. Counterstained with Hematoxylin.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myoglobin antibody [EP3081Y] (ab77232)


Purified ab77232 staining Myoglobin in Rat cardiac muscle tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with paraformaldehyde; antigen retrieval was by heat mediation in a Tris/EDTA buffer (**ab93684** Tris/EDTA buffer, pH 9.0). Samples were incubated with primary antibody (1/1000). A HRP-conjugated Goat Anti-Rabbit IgG H&L (ready to use) was used as the secondary antibody. Counterstained with Hematoxylin.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myoglobin antibody [EP3081Y] (ab77232)

Unpurified ab77232 showing positive staining in human Skeletal muscle tissue.


Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol.



Western blot - Anti-Myoglobin antibody [EP3081Y]
(ab77232)

This image is courtesy of an anonymous Abreview

Western blot analysis of Pig skeletal muscle whole tissue lysate (10µg/lane) labelling Myoglobin with unpurified ab77232 at 1/5000 in 1x PBS + 0.1% Tween-20 for 17 hours at 4°C. A IRDye® 680-conjugated Donkey anti-rabbit polyclonal (1/40000) was used as the secondary antibody.



Western blot - Anti-Myoglobin antibody [EP3081Y]
(ab77232)

All lanes : Anti-Myoglobin antibody [EP3081Y] (ab77232) at 1/5000 dilution (Unpurified)

Lane 1 : Human heart muscle lysate

Lane 2 : Mouse heart muscle lysate

Lane 3 : Rat heart muscle lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP labelled goat anti-rabbit at 1/2000 dilution

Predicted band size: 17 kDa


Observed band size: 17 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myoglobin antibody [EP3081Y] (ab77232)

Unpurified ab77232, at a 1/250 dilution, staining Myoglobin in paraffin embedded human heart muscle by Immunohistochemistry.


Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myoglobin antibody [EP3081Y] (ab77232)

Unpurified ab77232 showing negative staining in Normal human brain tissue.


Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myoglobin antibody [EP3081Y] (ab77232)

Unpurified ab77232 showing negative staining in Normal human kidney tissue.

Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Myoglobin antibody [EP3081Y] (ab77232)

Unpurified ab77232 showing negative staining in Normal human liver tissue.

Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol.

OI-RD Scanning - Anti-Myoglobin antibody
[EP3081Y] (ab77232)


Produced using unpurified ab77232.

Equilibrium disassociation constant (K_D)

Learn more about K_D

[Click here to learn more about \$K_D\$](#)

Why choose a recombinant antibody?



- Research with confidence**
Consistent and reproducible results
- Long-term and scalable supply**
Recombinant technology
- Success from the first experiment**
Confirmed specificity
- Ethical standards compliant**
Animal-free production

Anti-Myoglobin antibody [EP3081Y] (ab77232)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors