

# Anti-N6-methyladenosine (m6A) antibody ab151230

★★★★★ [3 Abreviews](#) [92 References](#) [画像数 1](#)

## 製品の概要

製品名	Anti-N6-methyladenosine (m6A) antibody
製品の詳細	Rabbit polyclonal to N6-methyladenosine (m6A)
由来種	Rabbit
アプリケーション	<b>適用あり:</b> Nucleotide Array
種交差性	<b>交差種:</b> Species independent
免疫原	Chemical/ Small Molecule corresponding to N6-methyladenosine (m6A) conjugated to keyhole limpet haemocyanin.
ポジティブ・コントロール	ab151230 was tested on a oligonucleotide array against RNA oligomers; N6-methyladenosine (m6A) and unmodified adenosine. 2Ome-RNA oligomers carry a methyl group at the 2'-OH residue of the ribose molecule, making them resistant to most ribonucleases.
特記事項	<p>N6-methyladenosine (m6A) is a post-transcriptional modification of RNA. m6A modification has been identified in all classes of RNA (rRNA, tRNA and mRNA) and is catalysed by an evolutionary conserved multi-subunit enzyme, methyltransferase like 3 (METTL3). Cellular and viral RNA has been known to be methylated for decades. Recent studies have found that mRNA is predominately m6A modified at stop codons and long internal exons, which are conserved between mouse and human. Emerging evidence suggests m6A plays an important role in regulating gene expression, alternative splicing patterns, downstream signalling (p53) as well as apoptosis.</p> <p>The regulation of m6A modifications in mRNA has been linked to disease, where fat mass and obesity-associated (FTO) has been has been reported to be a obesity risk gene. FTO is a m6A demethylase and polymorphisms that result in increased FTO expression are associated with increased body mass and risk of obesity.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at +4°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Constituent: PBS  PBS prepared using DEPC-treated water. The use of RNase inhibitor is recommended.
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab151230の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Nucleotide Array		Use a concentration of 2 µg/ml.

## ターゲット情報

**関連性**      N6-Methyladenosine (m6A) is an abundant modification found in mRNA, tRNA, snRNA, as well as long non-coding RNA, in all species. RNA adenosine methylation is catalyzed by a multicomponent complex composed of METTL3/MT-A70, METTL14, and WTAP in mammals. METTL3 & METTL14 are responsible for the methyltransferase activity of the complex, and WTAP mediates substrate recruitment.

## 画像



All batches of ab151230 are tested in Nucleotide Array against N6-methyladenosine (m6A) and unmodified adenosine. Six dilutions of each oligomer are printed on to the Array in triplicate and results are averaged before being plotted on to a graph. Results show strong binding to N6-methyladenosine, indicating that this antibody specifically recognises adenosine methylated at position N6.

2Ome-m6A - N6-methyladenosine

2Ome-m6A\_unmod - unmodified adenosine

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors