

# Anti-NCKAP1L antibody ab272631

画像数 5

### 製品の概要

<b>製品名</b>	Anti-NCKAP1L antibody
<b>製品の詳細</b>	Rabbit polyclonal to NCKAP1L
<b>由来種</b>	Rabbit
<b>アプリケーション</b>	<b>適用あり:</b> IHC-P, WB
<b>種交差性</b>	<b>交差種:</b> Human
<b>免疫原</b>	Recombinant fragment corresponding to Human NCKAP1L aa 1-100. Database link: <a href="#">P55160</a>
	<a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
<b>ポジティブ・コントロール</b>	WB: RT4, U-251MG sp, Human plasma (IgG/HSA depleted), Human liver tissue and tonsil tissue. IHC-P: Human lymph tissue, small intestine tissue, liver tissue, lung tissue.
<b>特記事項</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

<b>製品の状態</b>	Liquid
<b>保存方法</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>バッファー</b>	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), PBS
<b>精製度</b>	Immunogen affinity purified
<b>ポリ/モノ</b>	ポリクローナル
<b>アイソタイプ</b>	IgG

## アプリケーション


**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab272631の使用に適用されず  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use a concentration of 0.04 - 0.4 µg/ml.

## ターゲット情報

組織特異性	Expressed only in cells of hematopoietic origin.
配列類似性	Belongs to the HEM-1/HEM-2 family.
細胞内局在	Cell membrane.

## 画像



**Lanes 2-6 :** Anti-NCKAP1L antibody (ab272631) at 0.4 µg/ml

**Lane 1 :** Molecular weight ladder


**Lane 2 :** RT4

**Lane 3 :** U-251MG sp

**Lane 4 :** human plasma (IgG/HSA depleted)


**Lane 5 :** Human liver tissue

**Lane 6 :** Human tonsil tissue




Paraffin-embedded formalin-fixed Human lymph tissue stained for NCKAP1L using ab272631 at 1/500 in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NCKAP1L antibody (ab272631)




Paraffin-embedded formalin-fixed Human small intestine tissue stained for NCKAP1L using ab272631 at 1/500 in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NCKAP1L antibody (ab272631)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NCKAP1L antibody (ab272631)

Paraffin-embedded formalin-fixed Human liver tissue stained for NCKAP1L using ab272631 at 1/500 in immunohistochemical analysis.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-NCKAP1L antibody (ab272631)

Paraffin-embedded formalin-fixed Human lung tissue stained for NCKAP1L using ab272631 at 1/500 in immunohistochemical analysis.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors