

Anti-NKX2.8 antibody [NKX28/3233R] - BSA and Azide free ab270292

リコンビナント

画像数 2

製品の概要

製品名	Anti-NKX2.8 antibody [NKX28/3233R] - BSA and Azide free
製品の詳細	Rabbit monoclonal [NKX28/3233R] to NKX2.8 - BSA and Azide free
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment within Human NKX2.8 aa 10-123. The exact sequence is proprietary. Database link: O15522
ポジティブ・コントロール	WB: PC3 and MCF-7 cell lysate.
特記事項	ab270292 is a carrier free version of ab270261 . Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency. This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications. Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Constituent: PBS
キャリア・フリー	はい
精製度	Protein A/G purified

特記事項 (精製)	Purified from bioreactor concentrate.
ポリ/モノ	モノクローナル
クローン名	NKX28/3233R
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab270292の使用に適用されず**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 26 kDa.

ターゲット情報

配列類似性	Belongs to the NK-2 homeobox family. Contains 1 homeobox DNA-binding domain.
細胞内局在	Nucleus.

画像



All lanes : Anti-NKX2.8 antibody [NKX28/3233R] (**ab270261**) at 2 µg/ml

Lane 1 : PC3 (Human prostate adenocarcinoma cell line) cell lysate

Lane 2 : MCF-7 (Human breast adenocarcinoma cell line) cell lysate

Predicted band size: 26 kDa

This data was developed using the same antibody clone in a different buffer formulation containing PBS, BSA and sodium azide (**ab270261**).

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-NKX2.8 antibody [NKX28/3233R] - BSA and Azide free (ab270292)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors