

Anti-Nucleoporin p62/NUP62 antibody ab96134

★★★★★ 4 Abreviews 9 References 画像数 6

医薬用外劇物

製品の概要

製品名	Anti-Nucleoporin p62/NUP62 antibody
製品の詳細	Rabbit polyclonal to Nucleoporin p62/NUP62
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB, ICC/IF
種交差性	交差種: Mouse, Rat, Human 交差が予測される動物種: Cow, Zebrafish, Rhesus monkey 
免疫原	Recombinant fragment corresponding to Human Nucleoporin p62/NUP62 aa 294-502. Database link: NP_057637
ポジティブ・コントロール	WB: A431, H1299, PC-12, Rat2, Neuro-2a, GL261, C8D30, NIH/3T3, BCL-1, RAW 264.7 and C2C12 whole cell lysate. ICC/IF: A431 cells. IHC-P: Mouse fore brain tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab96134の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/1000.
WB	★★★★★ (2)	1/500 - 1/3000. Predicted molecular weight: 53 kDa.
ICC/IF	★★★★★ (1)	1/100 - 1/200.

ターゲット情報


関連性

The nuclear pore complex is a structure that extends across the nuclear envelope and regulates the flow of macromolecules between the cytoplasm and the nucleus. Nucleoporins are the main components of the nuclear pore complex in eukaryotic cells. Nup 62 is localized to the nuclear pore central plug. This protein associates with the importin alpha/beta complex which is involved in the import of proteins containing nuclear localization signals. There are multiple transcript variants of this gene, however, they encode a single protein isoform. This protein undergoes post-translational modifications. It contains about 10 N-acetylglucosamine side chain sites.

細胞内局在

Nucleus; nuclear pore complex. Cytoplasm; cytoskeleton; spindle pole. Note: Central region of the nuclear pore, within the transporter. During mitotic cell division, it associates with the poles of the mitotic spindle.

画像



All lanes : Anti-Nucleoporin p62/NUP62 antibody (ab96134) at 1/1000 dilution

Lane 1 : Neuro-2a (Mouse neuroblastoma cell line) whole cell lysate

Lane 2 : GL261 whole cell lysate

Lane 3 : C8D30 whole cell lysate

Lane 4 : NIH/3T3 (mouse embryo fibroblast cell line) whole cell lysate

Lane 5 : BCL-1 (mouse B lymphocyte lymphoma cell line) whole cell lysate

Lane 6 : RAW 264.7 (Mouse macrophage cell line transformed with Abelson murine leukemia virus) whole cell lysate

Lane 7 : C2C12 (Mouse myoblast cell line) whole cell lysate


Lysates/proteins at 30 µg per lane.

Secondary

All lanes : HRP-conjugated anti-rabbit IgG


Predicted band size: 53 kDa

10% SDS-PAGE.



Immunocytochemistry/ Immunofluorescence - Anti-Nucleoporin p62/NUP62 antibody (ab96134)

ab96134, at a 1/200 dilution, staining Nucleoporin p62/NUP62 in paraformaldehyde fixed A431 by Immunofluorescence analysis. Lower image: merged with DNA probe.



Western blot - Anti-Nucleoporin p62/NUP62 antibody (ab96134)

All lanes : Anti-Nucleoporin p62/NUP62 antibody (ab96134) at 1/1000 dilution

Lane 1 : PC-12 (Rat adrenal gland pheochromocytoma cell line) whole cell lysate

Lane 2 : Rat2 (Rat fibroblast cell line) whole cell lysate


Lysates/proteins at 30 µg per lane.

Secondary

All lanes : HRP-conjugated anti-rabbit IgG


Predicted band size: 53 kDa

10% SDS-PAGE.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Nucleoporin p62/NUP62 antibody (ab96134)

Paraffin-embedded mouse fore brain tissue stained for Nucleoporin p62/NUP62 using ab96134 at 1/500 dilution in immunohistochemical analysis.



Western blot - Anti-Nucleoporin p62/NUP62 antibody (ab96134)

All lanes : Anti-Nucleoporin p62/NUP62 antibody (ab96134) at 1/1000 dilution


Lane 1 : A431 whole cell lysate

Lane 2 : H1299 whole cell lysate

Lysates/proteins at 30 μ g per lane.

Predicted band size: 53 kDa

10% SDS PAGE



Immunofluorescence analysis of A431 cells, staining Nucleoporin p62/NUP62 with ab96134 at 1/200 dilution, demonstrating punctate nuclear staining.

Immunocytochemistry/ Immunofluorescence - Anti-Nucleoporin p62/NUP62 antibody (ab96134)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors