

Anti-PADI4 / PAD4 antibody [OTI4H5] ab128086

★★★★★ [1 Abreviews](#) [31 References](#) [画像数 13](#)

製品の概要

製品名	Anti-PADI4 / PAD4 antibody [OTI4H5]
製品の詳細	Mouse monoclonal [OTI4H5] to PADI4 / PAD4
由来種	Mouse
アプリケーション	適用あり: WB, IHC-P, ICC/IF, Flow Cyt (Intra)
種交差性	交差種: Human, African green monkey
免疫原	Recombinant full length protein, corresponding to amino acids 1-663 of Human PADI4/ PAD4 produced in HEK293T cells (NP_036519).
ポジティブ・コントロール	HEK293T cells transfected with PADI4 / PAD4 , Human breast, Carcinoma of Human lung, Adenocarcinoma of Human ovary, Human pancreas, Adenocarcinoma of Human endometrium, Human prostate, Carcinoma of Human prostate, Carcinoma of Human bladder and Human lymphoma tissues, COS7 cells transiently transfected by pCMV6-ENTRY PADI4 / PAD4 , Hela and Jurkat cells
特記事項	<p>The clone number has been updated from 4H5 to OTI4H5, both clone numbers name the same antibody clone.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	<p>pH: 7.30</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 48% PBS, 1% BSA, 50% Glycerol</p>
精製度	Protein A purified
特記事項 (精製)	Purified from TCS

ポリ/モノ	モノクローナル
クローン名	OT4H5
アイソタイプ	IgG2a

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab128086の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/2000. Predicted molecular weight: 74 kDa.
IHC-P	★★★★★ (1)	1/150.
ICC/IF		1/100.
Flow Cyt (Intra)		1/100. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

ターゲット情報

機能	Catalyzes the citrullination/deimination of arginine residues of proteins. Citrullinates histone H3 at 'Arg-8' and/or 'Arg-17' and histone H4 at 'Arg-3', which prevents their methylation by CARM1 and HRMT1L2/PRMT1 and represses transcription. Citrullinates EP300/P300 at 'Arg-2142', which favors its interaction with NCOA2/GRIP1.
組織特異性	Expressed in eosinophils and neutrophils, not expressed in peripheral monocytes or lymphocytes.
関連疾患	Genetic variations in PADI4 are a cause of susceptibility to rheumatoid arthritis (RA) [MIM:180300]. It is a systemic inflammatory disease with autoimmune features and a complex genetic component. It primarily affects the joints and is characterized by inflammatory changes in the synovial membranes and articular structures, widespread fibrinoid degeneration of the collagen fibers in mesenchymal tissues, and by atrophy and rarefaction of bony structures. Note=Could have an important role in the pathogenesis of rheumatoid arthritis by increasing citrullination of proteins in rheumatoid arthritis synovial tissues, leading, in a cytokine-rich milieu, to a break in tolerance to citrullinated peptides processed and presented in the appropriate HLA context.
配列類似性	Belongs to the protein arginine deiminase family.
細胞内局在	Cytoplasm. Nucleus. Cytoplasmic granule. Cytoplasmic granules of eosinophils and neutrophils.

画像



Western blot - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

All lanes : Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086) at 1/2000 dilution


Lane 1 : HEK293T cells transfected with pCMV6-ENTRY control cDNA

Lane 2 : HEK293T cells transfected with pCMV6-ENTRY PAD4 / PAD4 cDNA

Lysates/proteins at 5 µg per lane.


Predicted band size: 74 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC206501)




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Human breast tissue by Immunohistochemistry.




Immunocytochemistry/ Immunofluorescence - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/100, staining COS7 cells transiently transfected by pCMV6-ENTRY PAD4 / PAD4 by Immunofluorescence.




Flow Cytometry (Intracellular) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/100, staining PAD4 / PAD4 in HeLa cells (Red), compared to a nonspecific negative control antibody (Blue), analysed by Flow cytometry (Intracellular).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded carcinoma of Human lung tissue by Immunohistochemistry.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Adenocarcinoma of Human ovary tissue by Immunohistochemistry.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Human pancreas tissue by Immunohistochemistry.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Adenocarcinoma of Human endometrium tissue by Immunohistochemistry.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Human prostate tissue by Immunohistochemistry.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Carcinoma of Human prostate tissue by Immunohistochemistry.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Carcinoma of Human bladder tissue by Immunohistochemistry.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PAD4 / PAD4 antibody [OT14H5] (ab128086)

ab128086, at a dilution of 1/150, staining PAD4 / PAD4 in paraffin-embedded Human lymphoma tissue by Immunohistochemistry.



ab128086, at a dilution of 1/100, staining PAD14 / PAD4 in Jurkat cells (Red), compared to a nonspecific negative control antibody (Blue), analysed by Flow cytometry (Intracellular).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors