

Anti-PBR antibody ab2305

面格粒 1

医薬用外劇物

製品の概要

製品名	Anti-PBR antibody
製品の詳細	Rabbit polyclonal to PBR
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Jurkat cell lysate. Recombinant PBR. Rat kidney tissue lysate.
特記事項	This product is manufactured by BioVision, an Abcam company and was previously called 3044 PBR Antibody. 3044-100 is the same size as the 100 µg size of ab2305.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 30% Glycerol (glycerin, glycerine), 0.5% BSA, 69% PBS
精製度	Affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab2305の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 - 4 µg/ml.

ターゲット情報

機能	Responsible for the manifestation of peripheral-type benzodiazepine recognition sites and is most likely to comprise binding domains for benzodiazepines and isoquinoline carboxamides. May play a role in the transport of porphyrins and heme. Plays a role in the transport of cholesterol across mitochondrial membranes in steroidogenic cells.
組織特異性	Found in many tissue types. Expressed at the highest levels under normal conditions in tissues that synthesize steroids.
配列類似性	Belongs to the TspO/BZRP family.
細胞内局在	Mitochondrion membrane.

画像



All lanes : Anti-PBR antibody (ab2305) at 4 µg/ml

Lane 1 : Jurkat (Human T cell leukemia cell line from peripheral blood) cell lysate

Lane 2 : Recombinant PBR

Lane 3 : Rat kidney tissue lysate

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors