

Anti-PDE1C/PDE1 antibody ab14602

7 References [画像数 1](#)

製品の概要

製品名	Anti-PDE1C/PDE1 antibody
製品の詳細	Rabbit polyclonal to PDE1C/PDE1
由来種	Rabbit
特異性	The antibody does not cross react with PDE1A, PDE1B or other PDEs.
アプリケーション	適用あり: IP, IHC-P, WB
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab14602の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration. PubMed: 19924104
WB		1/250 - 1/500. Predicted molecular weight: 70 kDa. This antibody labels multiple bands in the 70-87 kDa range all believed to be various PDE1C variants.

ターゲット情報

機能	Cyclic nucleotide phosphodiesterase with a dual-specificity for the second messengers cAMP and cGMP, which are key regulators of many important physiological processes. Has a high affinity for both cAMP and cGMP.
組織特異性	Expressed in several tissues, including brain and heart.
配列類似性	Belongs to the cyclic nucleotide phosphodiesterase family. PDE1 subfamily.

画像



All lanes : Anti-PDE1C/PDE1 antibody (ab14602) at 1/500 dilution

Lane 1 : Brain lysate

Lane 2 : Heart lysate

Lane 3 : Lung lysate

Lane 4 : Testis lysate

Lane 5 : Kidney lysate

Predicted band size: 70 kDa

Western Blot analysis of PDE1C/PDE1 with ab14602 (1/500). Rat tissues: 1- Brain; 2 -Heart; 3 - Lung; 4 - Testis; 5 - Kidney.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors