

### PE Anti-CD30 antibody [2SH12-5F] ab25187

画像数 1

#### 製品の概要

製品名	PE Anti-CD30 antibody [2SH12-5F]
製品の詳細	PE Armenian hamster monoclonal [2SH12-5F] to CD30
由来種	Armenian hamster
標識	PE. Ex: 488nm, Em: 575nm
アプリケーション	<b>適用あり:</b> Flow Cyt
種交差性	<b>交差種:</b> Mouse
免疫原	Recombinant CD30 extracellular domain-mouse IgG1 fusion protein (Mouse)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	<p>pH: 7.3</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 16% Sucrose</p> <p>Stabilising agent.</p>
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	2SH12-5F
アイソタイプ	IgG
軽鎖の種類	kappa

## アプリケーション

### The Abpromise guarantee

**Abpromise保証は、**次のテスト済みアプリケーションにおけるab25187の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 0.2µg for 10 <sup>6</sup> cells.

## ターゲット情報

### 機能

Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B.

### 配列類似性

Contains 6 TNFR-Cys repeats.


### 翻訳後修飾

Phosphorylated on serine and tyrosine residues.

### 細胞内局在

Cytoplasm and Cell membrane.

## 画像



Flow cytometry analysis staining CD30 in NIH/3T3 and mCD30-Ig transfected 3T3 cells using ab25187 at a dilution of 0.1 µg/10<sup>6</sup> cells. Hamster IgG-PE was used as an isotype control.

Flow Cytometry - PE Anti-CD30 antibody [2SH12-5F] (ab25187)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors