

PE Mouse IgG2b [PLPV219] - Isotype Control ab91532

4 References

製品の概要

| | |
|----------|--|
| 製品名 | PE Mouse IgG2b [PLPV219] - Isotype Control |
| 標識 | PE. Ex: 488nm, Em: 575nm |
| 特異性 | ab91532 was produced against a synthetic hapten, which is normally not present in Humans or animals. Therefore, ab91532 does not react with any antigen of Human or animal cells. |
| アプリケーション | 適用あり: Flow Cyt |
| 免疫原 | Vegetal protein |
| 特記事項 | <p>25µg is sufficient for 50 tests</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. Store In the Dark. |
| バッファー | <p>pH: 7.20</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituents: PBS, 1% BSA</p> |
| 精製度 | Protein G purified |
| ポリ/モノ | モノクローナル |
| クローン名 | PLPV219 |
| アイソタイプ | IgG2b |

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab91532の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| Flow Cyt | | |

追加情報

Flow Cyt: Use 20µl for 10⁶ cells.
ab91532 is a negative isotype control

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors