

PerCP Anti-CD80 antibody [MEM-233] ab82410

1 References [画像数 1](#)

製品の概要

製品名	PerCP Anti-CD80 antibody [MEM-233]
製品の詳細	PerCP Mouse monoclonal [MEM-233] to CD80
由来種	Mouse
標識	PerCP. Ex: 482nm, Em: 675nm
特異性	ab82410 reacts with CD80 expressed on professional antigen presenting cells, such as dendritic cells, macrophages or activated B lymphocytes.
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	Extracellular domain of human CD80 fused to human IgG1(Fc)
ポジティブ・コントロール	RAJI Human Burkitt lymphoma cell line.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.4</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituent: 0.2% BSA</p> <p>0.2% high-grade protease free BSA</p>
精製度	Size exclusion
特記事項(精製)	Conjugate made from highly purified antibody (>95% by SDS-PAGE).
ポリ/モノ	モノクローナル

クローン名	MEM-233
アイソタイプ	IgG1

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab82410の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		<p>Use 10µl for 10⁶ cells. or 100 µl of whole blood.</p> <p>ab118652 Mouse monoclonal IgG1 is suitable for use as an</p>

ターゲット情報

機能	Involved in the costimulatory signal essential for T-lymphocyte activation. T-cell proliferation and cytokine production is induced by the binding of CD28 or CTLA-4 to this receptor.
組織特異性	Expressed on activated B-cells, macrophages and dendritic cells.
配列類似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
細胞内局在	Membrane.

画像



Flow Cytometry analysis of P815 cells labeling CD80 with Anti-CD80 antibody [MEM-233] (PerCP) (ab82410).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors