

PerCP/Cy5.5® Anti-CD3 antibody [UCHT1], prediluted ab157317

4 References [画像数 1](#)

製品の概要

製品名	PerCP/Cy5.5® Anti-CD3 antibody [UCHT1], prediluted
製品の詳細	PerCP/Cy5.5® Mouse monoclonal [UCHT1] to CD3, prediluted
由来種	Mouse
標識	PerCP/Cy5.5®. Ex: 482nm, Em: 690nm
特異性	ab157317 reacts with the epsilon chain of the CD3 complex.
アプリケーション	適用あり: Flow Cyt, Flow Cyt (Intra)
種交差性	交差種: Human, Non human primates
免疫原	Tissue, cells or virus corresponding to Human CD3. Human thymocytes followed by Sezary T cells.
ポジティブ・コントロール	Flow Cytometry: human peripheral blood cells. Flow Cyt (intra): Human blood cells.
特記事項	<p>Purified antibody is conjugated with tandem dye PerCP/Cy5.5® under optimum conditions. The conjugate is purified by size-exclusion chromatography and adjusted for direct use. No reconstitution is necessary.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	<p>pH: 7.4</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituents: 99% PBS, 0.2% BSA</p>
精製度	Size exclusion
ポリ/モノ	モノクローナル


クローン名	UCHT1
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab157317の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 4µl for 10 ⁶ cells. (or 100 µl of whole blood). ab157226 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
Flow Cyt (Intra)		Use 4µl for 10 ⁶ cells. (or 100 µl of whole blood). ab157226 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

画像



Analysis by flow cytometry of ab157317 staining CD3 in human peripheral blood cells.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors