

### Anti-Perilipin 2 antibody [ADFP/1494] ab219301

#### 画像数 1

#### 製品の概要

---

製品名	Anti-Perilipin 2 antibody [ADFP/1494]
製品の詳細	Mouse monoclonal [ADFP/1494] to Perilipin-2
由来種	Mouse
アプリケーション	<b>適用あり:</b> IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment within Human Perilipin-2 aa 200-400. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please <b>contact</b> our Scientific Support team to discuss your requirements. Database link: <a href="#">Q99541</a>
ポジティブ・コントロール	Human adrenal tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.05% Sodium azide Constituents: 99% PBS, 0.05% BSA
精製度	Protein A/G purified
特記事項(精製)	ab219301 was purified from Bioreactor Concentrate by Protein A/G.
ポリ/モノ	モノクローナル
クローン名	ADFP/1494

アイソタイプ	IgG2b
軽鎖の種類	kappa

## アプリケーション


**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab219301の使用に適用されず  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 2 - 4 µg/ml. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol. (Primary incubation for 30 minutes at RT).

## ターゲット情報

機能	May be involved in development and maintenance of adipose tissue.
組織特異性	Milk lipid globules.
配列類似性	Belongs to the perilipin family.
翻訳後修飾	Acylated; primarily with C14, C16 and C18 fatty acids.
細胞内局在	Membrane.

## 画像



Immunohistochemical analysis of formalin-fixed, paraffin-embedded human adrenal tissue labeling Perilipin-2 with ab219301 at 4 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Perilipin 2 antibody [ADFP/1494] (ab219301)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors