

### Anti-Peripherin antibody ab39374

★★★★★ [7 Abreviews](#) [19 References](#) [画像数 1](#)

#### 製品の概要

製品名	Anti-Peripherin antibody
製品の詳細	Chicken polyclonal to Peripherin
由来種	Chicken
アプリケーション	<b>適用あり:</b> ICC/IF
種交差性	<b>交差種:</b> Rat
免疫原	Recombinant protein expressed in bacteria.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: PBS, 1% BSA</p> <p>(10mM isotonic 0.9% w/v), pH 7.2</p>
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

#### アプリケーション

## The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab39374の使用に適用されます**

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.

## ターゲット情報

**機能** Class-III neuronal intermediate filament protein.

**配列類似性** Belongs to the intermediate filament family.

## 画像



Ab39374 staining human peripherin in neonatal rat cortical neurons cultured for 14-21 days by ICC/IF.

Immunocytochemistry/ Immunofluorescence - Anti-Peripherin antibody (ab39374)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors