

Anti-PKC antibody [EPR16898] ab179522

リコンビナント **RabMAb**

8 References [画像数 12](#)

製品の概要

製品名	Anti-PKC antibody [EPR16898]
製品の詳細	Rabbit monoclonal [EPR16898] to PKC
由来種	Rabbit
アプリケーション	適用あり: Flow Cyt (Intra), WB, ICC/IF, IP
種交差性	交差種: Mouse, Rat, Human
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: HeLa, 293T, SH-SY5Y, K562, A431, L6, Neuro-2a and L929 whole cell lysates; Human fetal brain, fetal heart and fetal kidney lysates; Rat brain whole cell lysate. ICC/IF: K562 cells. Flow Cyt (intra): NIH/3T3 cells. IP: Rat brain extract.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR16898
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab179522の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/360. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
WB		1/1000. Detects a band of approximately 77 kDa (predicted molecular weight: 77 kDa).
ICC/IF		1/200.
IP		1/100.

ターゲット情報

機能

Calcium-activated, phospholipid- and diacylglycerol (DAG)-dependent serine/threonine-protein kinase that is involved in positive and negative regulation of cell proliferation, apoptosis, differentiation, migration and adhesion, tumorigenesis, cardiac hypertrophy, angiogenesis, platelet function and inflammation, by directly phosphorylating targets such as RAF1, BCL2, CSPG4, TNNT2/CTNT, or activating signaling cascade involving MAPK1/3 (ERK1/2) and RAP1GAP. Involved in cell proliferation and cell growth arrest by positive and negative regulation of the cell cycle. Can promote cell growth by phosphorylating and activating RAF1, which mediates the activation of the MAPK/ERK signaling cascade, and/or by up-regulating CDKN1A, which facilitates active cyclin-dependent kinase (CDK) complex formation in glioma cells. In intestinal cells stimulated by the phorbol ester PMA, can trigger a cell cycle arrest program which is associated with the accumulation of the hyper-phosphorylated growth-suppressive form of RB1 and induction of the CDK inhibitors CDKN1A and CDKN1B. Exhibits anti-apoptotic function in glioma cells and protects them from apoptosis by suppressing the p53/TP53-mediated activation of IGFBP3, and in leukemia cells mediates anti-apoptotic action by phosphorylating BCL2. During macrophage differentiation induced by macrophage colony-stimulating factor (CSF1), is translocated to the nucleus and is associated with macrophage development. After wounding, translocates from focal contacts to lamellipodia and participates in the modulation of desmosomal adhesion. Plays a role in cell motility by phosphorylating CSPG4, which induces association of CSPG4 with extensive lamellipodia at the cell periphery and polarization of the cell accompanied by increases in cell motility. Is highly expressed in a number of cancer cells where it can act as a tumor promoter and is implicated in malignant phenotypes of several tumors such as gliomas and breast cancers. Negatively regulates myocardial contractility and positively regulates angiogenesis, platelet aggregation and thrombus formation in arteries. Mediates hypertrophic growth of neonatal cardiomyocytes, in part through a MAPK1/3 (ERK1/2)-dependent signaling pathway, and upon PMA treatment, is required to induce cardiomyocyte hypertrophy up to heart failure and death, by increasing protein synthesis, protein-DNA ratio and cell surface area. Regulates cardiomyocyte function by phosphorylating cardiac troponin T (TNNT2/CTNT), which induces significant reduction in actomyosin ATPase activity, myofilament calcium sensitivity and myocardial contractility. In angiogenesis, is required for full endothelial cell migration, adhesion to vitronectin (VTN), and vascular endothelial growth factor A (VEGFA)-dependent regulation of kinase activation and vascular tube formation. Involved in the stabilization of VEGFA mRNA at

post-transcriptional level and mediates VEGFA-induced cell proliferation. In the regulation of calcium-induced platelet aggregation, mediates signals from the CD36/GP4 receptor for granule release, and activates the integrin heterodimer ITGA2B-ITGB3 through the RAP1GAP pathway for adhesion. During response to lipopolysaccharides (LPS), may regulate selective LPS-induced macrophage functions involved in host defense and inflammation. But in some inflammatory responses, may negatively regulate NF-kappa-B-induced genes, through IL1A-dependent induction of NF-kappa-B inhibitor alpha (NFKBIA/IKBA). Upon stimulation with 12-O-tetradecanoylphorbol-13-acetate (TPA), phosphorylates EIF4G1, which modulates EIF4G1 binding to MKNK1 and may be involved in the regulation of EIF4E phosphorylation. Phosphorylates KIT, leading to inhibition of KIT activity. Phosphorylates ATF2 which promotes cooperation between ATF2 and JUN, activating transcription.


配列類似性

Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. PKC subfamily. Contains 1 AGC-kinase C-terminal domain.
Contains 1 C2 domain.
Contains 2 phorbol-ester/DAG-type zinc fingers.
Contains 1 protein kinase domain.

細胞内局在

Cytoplasm. Cell membrane. Mitochondrion membrane. Nucleus.

画像



All lanes : Anti-PKC antibody [EPR16898] (ab179522) at 1/2000 dilution

Lane 1 : Active human PKC alpha full length protein

Lane 2 : Active human PKC beta 1 full length protein

Lane 3 : Active human PKC beta 2 full length protein

Lane 4 : Active human PKC gamma full length protein

Lane 5 : Active human PKC delta full length protein

Lane 6 : Active human PKC eta full length protein

Lane 7 : Active human PKC epsilon full length protein

Lane 8 : Active human PKC theta full length protein

Lane 9 : Active human PKC mu full length protein

Lysates/proteins at 0.02 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#))


Predicted band size: 77 kDa

Exposure time: 2 seconds

Active human PKC alpha full length protein (Catalog#[ab55672](#))

contains aa1-672 with GST-tag; Active human PKC beta 1 full length protein(Catalog#[ab60840](#)) contains aa1-671 with GST-tag; Active human PKC beta 2 full length protein(Catalog#[ab60841](#)) contains aa1-673 with GST-tag; Active human PKC gamma full length protein(Catalog#[ab60842](#)) contains aa1-677 with GST-tag; Active human PKC delta full length protein(Catalog#[ab60844](#)) contains aa1-676 with GST-tag; Active human PKC eta full length protein(Catalog#[ab60849](#)) contains aa1-683 with GST-tag; Active human PKC epsilon full length protein(Catalog#[ab60847](#)) contains aa1-737 with GST-tag; Active human PKC theta full length protein(Catalog#[ab56641](#)) contains aa1-706 with GST-tag; Active human PKC mu full length protein(Catalog#[ab60873](#)) contains aa1-912 with GST-tag.


Blocking/Diluting buffer 5% NFDM/TBST



Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized K562 (Human chronic myelogenous leukemia cells from bone marrow) cells labeling PKC alpha + beta + gamma with ab179522 at 1/200 dilution, followed by Goat **anti-rabbit Alexa Fluor® 488** (IgG) ([ab150077](#)) secondary antibody at 1/400 dilution (green). Confocal image shows cytoplasmic staining on K562 cell line. The nuclear counterstain is DAPI (blue). Tubulin is detected with [ab7291](#) (anti-Tubulin mouse mAb) at 1/500 dilution and [ab150120](#) (goat **anti-mouse AlexaFluor®594** secondary antibody) at 1/500 dilution (red).

The negative controls are as follows:-

-ve control 1: - ab179522 at 1/200 dilution followed by [ab150120](#) (AlexaFluor®594 Goat anti-Mouse secondary) at 1/500 dilution.
 -ve control 2: - [ab7291](#) (anti-Tubulin mouse mAb) at 1/500 dilution followed by [ab150077](#) (Alexa Fluor®488 Goat Anti-Rabbit IgG H&L) at 1/400 dilution.



Western blot - Anti-PKC antibody [EPR16898]
(ab179522)

All lanes : Anti-PKC antibody [EPR16898] (ab179522) at 1/2000 dilution

Lane 1 : Active human PKC zeta full length protein

Lane 2 : Active human PKC iota full length protein

Lysates/proteins at 0.02 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/50000 dilution

Predicted band size: 77 kDa

Exposure time: 3 minutes


Active human PKC zeta full length protein (Catalog# **ab60848**) contains aa1-592 with GST-tag; Active human PKC iota full length protein (Catalog# **ab60850**) contains aa1-596 with GST-tag.

Blocking/Diluting buffer 5% NFDM/TBST

PKC alpha	PKC beta I	PKC beta II	PKC gamma	PKC delta	PKC epsilon	PKC eta	PKC mu	PKC theta	PKC iota	PKC zeta
+	+	+	+	-	+	+	-	+	-	-

Western blot - Anti-PKC antibody [EPR16898]
(ab179522)

Recombinant PKC isoforms were tested for reactivity by western blot, and are marked as positive (+) or negative (-).



Western blot - Anti-PKC antibody [EPR16898]
(ab179522)

All lanes : Anti-PKC antibody [EPR16898] (ab179522) at 1/10000 dilution

Lane 1 : HeLa (Human epithelial cells from cervix adenocarcinoma) whole cell lysate at 20 µg

Lane 2 : 293T (Human epithelial cells from embryonic kidney) whole cell lysate at 10 µg

Lane 3 : SH-SY5Y (Human neuroblastoma from bone marrow cells) whole cell lysate at 10 µg


Secondary

All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 77 kDa

Observed band size: 77 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.



All lanes : Anti-PKC antibody [EPR16898] (ab179522) at 1/1000 dilution

Lane 1 : K562 (Human chronic myelogenous leukemia cells from bone marrow) whole cell lysate at 20 µg

Lane 2 : A431 (Human epidermoid carcinoma) whole cell lysate at 10 µg


Secondary

All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 77 kDa

Observed band size: 77 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.



All lanes : Anti-PKC antibody [EPR16898] (ab179522) at 1/1000 dilution

Lane 1 : Human fetal brain tissue lysate

Lane 2 : Human fetal heart tissue lysate

Lane 3 : Human fetal kidney tissue lysate

Lysates/proteins at 10 µg per lane.


Secondary

All lanes : Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at 1/1000 dilution

Predicted band size: 77 kDa

Observed band size: 77 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.



Western blot - Anti-PKC antibody [EPR16898] (ab179522)

All lanes : Anti-PKC antibody [EPR16898] (ab179522) at 1/5000 dilution

Lane 1 : L6 (Rat skeletal muscle cell line) whole cell lysate

Lane 2 : Neuro-2a (Mouse neuroblastoma cells) whole cell lysate

Lane 3 : L-929 (Mouse connective tissue fibroblast cells) whole cell lysate

Lysates/proteins at 10 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 77 kDa

Observed band size: 77 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.



Western blot - Anti-PKC antibody [EPR16898] (ab179522)

Anti-PKC antibody [EPR16898] (ab179522) at 1/1000 dilution + Rat brain whole tissue lysate at 10 µg


Secondary

Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 77 kDa


Observed band size: 77 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.



Flow Cytometry (Intracellular) - Anti-PKC antibody
[EPR16898] (ab179522)

Intracellular flow cytometric analysis of 2% paraformaldehyde-fixed NIH/3T3 (Mouse embryo fibroblast cells) cells labeling PKC alpha + beta + gamma with ab179522 at 1/360 dilution (red) compared with a rabbit monoclonal IgG isotype control (black) and an unlabelled control (cells without incubation with primary antibody and secondary antibody; blue). Goat anti rabbit IgG (FITC) at 1/150 dilution was used as the secondary antibody.



Immunoprecipitation - Anti-PKC antibody
[EPR16898] (ab179522)

PKC alpha + beta + gamma was immunoprecipitated from 1mg of Rat brain extract with ab179522 at 1/100 dilution. Western blot was performed from the immunoprecipitate using ab179522 at 1/1000 dilution. Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated antibody was used as secondary antibody at 1/1000 dilution.

Blocking and dilution buffer and concentration: 5% NFDM/TBST.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-PKC antibody [EPR16898] (ab179522)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors