

# Anti-Poliovirus (type 1-3) antibody ab22450

## 3 References

### 製品の概要

---

製品名	Anti-Poliovirus (type 1-3) antibody
製品の詳細	Goat polyclonal to Poliovirus (type 1-3)
由来種	Goat
特異性	Poliovirus Types 1, 2 & 3
アプリケーション	<b>適用あり:</b> ICC/IF
種交差性	<b>交差種:</b> Enterovirus C
免疫原	Tissue, cells or virus corresponding to Poliovirus (type 1-3). Pool of purified Poliovirus (types 1-3). The strains used in the production of this antibody are Mahoney, MEF-1, and Saukett.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.1% Sodium azide Constituent: PBS
精製度	Proprietary Purification
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

---

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab22450の使用に適用されず**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent dilution.

## ターゲット情報

---

**関連性**      Poliovirus is a single stranded RNA virus and a member of the Picornavirus family. It has three strains: Poliovirus 1, 2 and 3.

**細胞内局在**      Whole virus

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors