

Anti-PVL LukS antibody ab190473

1 References [画像数 1](#)

製品の概要

製品名	Anti-PVL LukS antibody
製品の詳細	Rabbit polyclonal to PVL LukS
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Staphylococcus aureus
免疫原	Recombinant full length protein (His-tag) corresponding to PVL LukS. (Staphylococcus aureus PVL LukS purified from E.coli).
ポジティブ・コントロール	WB: Purified PVL Luk-S.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituent: 99% PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab190473の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 µg/ml. Predicted molecular weight: 35 kDa. If probing against purified protein, at least 10 ng of pure toxin can be detected with 0.5 ug/ml of the primary antibody combined with a suitable dilution of a secondary anti-rabbit antibody

画像



Anti-PVL LukS antibody (ab190473) at 0.5 µg/ml + Purified PVL Luk-S

Secondary

Goat anti rabbit-AP at 1/3000 dilution

Predicted band size: 35 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors