

Anti-RAP1GAP antibody [Y134] ab32373

リコンビナント **RabMAb**

★★★★★ **1 Abreviews** **11 References** 画像数 **12**

製品の概要

製品名	Anti-RAP1GAP antibody [Y134]
製品の詳細	Rabbit monoclonal [Y134] to RAP1GAP
由来種	Rabbit
特異性	ab32373 recognises Rap1 GTPase-activating protein 1
アプリケーション	適用あり: Flow Cyt (Intra), ICC/IF, WB, IHC-P, IP
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Jurkat, SH-SY5Y cell lysate, mouse and rat brain lysate IHC: Human breast cancer ICC/IF: HeLa cells IP: C6 cell lysate
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	Y134

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab32373の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/30 - 1/50. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
ICC/IF		1/100 - 1/250.
WB		1/10000 - 1/50000. Detects a band of approximately 95 kDa (predicted molecular weight: 73 kDa).
IHC-P	★★★★★ (1)	1/1000. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.
IP		1/20 - 1/60.

ターゲット情報

機能	GTPase activator for the nuclear Ras-related regulatory protein RAP-1A (KREV-1), converting it to the putatively inactive GDP-bound state.
組織特異性	Significant expression seen in the brain, kidney and pancreas. Abundant in the cerebral cortex and expressed at much lower levels in the spinal cord. Not detected in the lymphoid tissues.
配列類似性	Contains 1 GoLoco domain. Contains 1 Rap-GAP domain.
細胞内局在	Golgi apparatus membrane.

画像



All lanes : Anti-RAP1GAP antibody [Y134] (ab32373) at 1/10000 dilution (purified)

Lane 1 : Mouse brain lysate

Lane 2 : Rat brain lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution


Predicted band size: 73 kDa

Observed band size: 95 kDa

Blocking buffer: 5% NFDM/TBST

Dilution buffer: 5% NFDM/TBST

The molecular weight observed is consistent with what has been described in the literature (PMID: 9346962, PMID: 15632203, PMID: 30144784).



Anti-RAP1GAP antibody [Y134] (ab32373) at 1/10000 dilution (purified) + Jurkat cell lysate at 20 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution


Predicted band size: 73 kDa

Observed band size: 95 kDa

Blocking buffer: 5% NFDM/TBST

Dilution buffer: 5% NFDM/TBST

The molecular weight observed is consistent with what has been described in the literature (PMID: 9346962, PMID: 15632203, PMID: 30144784).



Western blot - Anti-RAP1GAP antibody [Y134]
(ab32373)

Anti-RAP1GAP antibody [Y134] (ab32373) at 1/10000 dilution
(purified) + SH-SY5Y cell lysate at 20 µg

Secondary

Anti-rabbit IgG (HRP), specific to the non-reduced form of IgG at
1/1000 dilution


Predicted band size: 73 kDa

Observed band size: 95 kDa

Blocking buffer: 5% NFDm/TBST

Dilution buffer: 5% NFDm/TBST

The molecular weight observed is consistent with what has been
described in the literature (PMID: 9346962, PMID: 15632203,
PMID: 30144784).




Western blot - Anti-RAP1GAP antibody [Y134]
(ab32373)

Anti-RAP1GAP antibody [Y134] (ab32373) at 1/50000 dilution
(unpurified) + Jurkat cell lysate

Predicted band size: 73 kDa


Observed band size: 95 kDa

The molecular weight observed is consistent with what has been
described in the literature (PMID: 9346962, PMID: 15632203,
PMID: 30144784).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-RAP1GAP antibody [Y134] (ab32373)

Immunohistochemical staining of paraffin embedded human cerebral cortex with purified ab32373 at a working dilution of 1/1000. The secondary antibody used is **ab97051**, a goat anti-rabbit IgG (H&L) at a dilution of 1/500. The sample is counter-stained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0. PBS was used instead of the primary antibody as the negative control, and is shown in the inset.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-RAP1GAP antibody [Y134] (ab32373)

Immunohistochemical staining of paraffin-embedded human breast cancer tissue using unpurified ab32373 at 1/100 dilution.




Immunocytochemistry/ Immunofluorescence - Anti-RAP1GAP antibody [Y134] (ab32373)

Immunofluorescence staining of HeLa cells with purified ab32373 at a working dilution of 1/250, counter-stained with DAPI. The secondary antibody was Alexa Fluor® 488 goat anti-rabbit ([ab150077](#)), used at a dilution of 1/1000. [ab7291](#), a mouse anti-tubulin antibody (1/1000), was used to stain tubulin along with [ab150120](#) (Alexa Fluor® 594 goat anti-mouse, 1/1000), shown in the top right hand panel. The cells were fixed in 4 % PFA and permeabilized using 0.1% Triton X 100. The negative controls are shown in bottom middle and right hand panels - for negative control 1, purified ab32373 was used at a dilution of 1/500 followed by an Alexa Fluor® 594 goat anti-mouse antibody ([ab150120](#)) at a dilution of 1/500. For negative control 2, [ab7291](#) (mouse anti-tubulin) was used at a dilution of 1/500 followed by an Alexa Fluor® 488 goat anti-rabbit antibody ([ab150077](#)) at a dilution of 1/400.




Immunocytochemistry/ Immunofluorescence - Anti-RAP1GAP antibody [Y134] (ab32373)

ICC/IF image of unpurified ab32373 stained HeLa cells. The cells were 4% formaldehyde fixed (10 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab32373, 1µg/ml) overnight at +4°C. The secondary antibody (green) was [ab96899](#), DyLight® 488 goat anti-rabbit IgG (H+L) used at a 1/250 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.




Flow Cytometry (Intracellular) - Anti-RAP1GAP antibody [Y134] (ab32373)

Overlay histogram showing Jurkat cells fixed in 4% PFA and stained with purified ab32373 at a dilution of 1/30 (red line). The secondary antibody used was Alexa Fluor® 488 goat anti-rabbit at a dilution of 1/500. Rabbit monoclonal IgG was used as an isotype control (black line) and cells incubated in the absence of both primary and secondary antibody were used as a negative control (blue line).



Flow Cytometry (Intracellular) - Anti-RAP1GAP antibody [Y134] (ab32373)

Overlay histogram showing SH-SY5Y cells stained with unpurified ab32373 (red line). The cells were fixed with 4% paraformaldehyde (10 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab32373, 1/100) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-rabbit IgG (H+L) ([ab96899](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was rabbit IgG (monoclonal) (1µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.



ab32373 (purified) at 1/20 immunoprecipitating RAP1GAP in 10 µg C6 cell lysate (Lanes 1 and 2, observed at 82-95 kDa). Lane 3 - Rabbit monoclonal IgG (**ab172730**). For western blotting, HRP Veriblot for IP Detection Reagent (**ab131366**) was used for detection (1/1000). Blocking buffer and concentration: 5% NFDm/TBST Dilution buffer and concentration: 5% NFDm/TBST

Why choose a recombinant antibody?

Research with confidence
Consistent and reproducible results

Long-term and scalable supply
Recombinant technology

Success from the first experiment
Confirmed specificity

Ethical standards compliant
Animal-free production

Anti-RAP1GAP antibody [Y134] (ab32373)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors