# abcam

### Product datasheet

## Anti-RFC5 antibody ab3619

2 References 画像数1

製品の概要

製品名 Anti-RFC5 antibody

製品の詳細 Goat polyclonal to RFC5

由来種 Goat

特異性 We have data to indicate that this antibody may not cross react with Mouse. However, this has not

been conclusively tested and expression levels may vary in certain cell lines/tissues.

アプリケーション **適用あり:** WB

種交差性 交差種: Human

交差が予測される動物種: Chimpanzee, Gorilla, Orangutan

免疫原 Synthetic peptide corresponding to Human RFC5.

ポジティブ・コントロール HeLa cell nuclear extract.

特記事項 The Life Science industry has been in the grips of a reproducibility crisis for a number of years.

> Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets

your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be

found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態 Liquid

保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

バッファー pH: 7

Preservative: 0.1% Sodium azide

Constituents: 0.021% PBS, 1.764% Sodium citrate, 1.815% Tris

Immunogen affinity purified 精製度

ポリ/モノ ポリクローナル

アイソタイプ lgG **The Abpromise guarantee Abpromise保証は、**次のテスト済みアプリケーションにおけるab3619の使用に適用されます アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/5000. Detects a band of approximately 37 kDa (predicted molecular weight: 36 kDa).

#### ターゲット情報


機能 The elongation of primed DNA templates by DNA polymerase delta and epsilon requires the

action of the accessory proteins proliferating cell nuclear antigen (PCNA) and activator 1.

**配列類似性** Belongs to the activator 1 small subunits family.

細胞内局在 Nucleus.

#### 画像



Sample: 50µg Hela cell nuclear extract

Anti-RFC5 (ab3619) used at 1:1000 dilution

Sample: 50µg Hela cell nuclear extract

Anti-RFC5 (ab3619) used at 1:1000 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- · Replacement or refund for products not performing as stated on the datasheet
- · Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- · We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <a href="https://www.abcam.co.jp/abpromise">https://www.abcam.co.jp/abpromise</a> or contact our technical team.

#### Terms and conditions

· Guarantee only valid for products bought direct from Abcam or one of our authorized distributors