


Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] ab181845

リコンビナント **RabMAb**

★★★★★ **1 Abreviews** **4 References** 画像数 5

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] |
| 製品の詳細 | Rabbit monoclonal [EPR7266(2)] to SCDGFB/PDGF-D |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB, IHC-P |
| 種交差性 | 交差種: Mouse, Human 交差が予測される動物種: Rat  |
| 免疫原 | Synthetic peptide. This information is proprietary to Abcam and/or its suppliers. |
| ポジティブ・コントロール | Human fetal kidney, fetal heart and small intestine lysate. PC-3 cell lysate. Human squamous cell lung carcinoma tissue. Mouse heart tissue. |
| 特記事項 | <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| バッファー | Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA |
| 精製度 | Protein A purified |
| ポリモノ | モノクローナル |
| クローン名 | EPR7266(2) |
| アイソタイプ | IgG |

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab181845の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | 1/1000 - 1/2000. Predicted molecular weight: 43 kDa. |
| IHC-P | ★★★★★ (1) | 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |

ターゲット情報

| | |
|-------|---|
| 機能 | Potent mitogen for cells of mesenchymal origin. Binding of this growth factor to its affinity receptor elicits a variety of cellular responses. It is released by platelets upon wounding and plays an important role in stimulating adjacent cells to grow and thereby heals the wound. Activated by proteolytic cleavage and this active form acts as a specific ligand for beta platelet-derived growth factor receptor. Induces macrophage recruitment, increased interstitial pressure, and blood vessel maturation during angiogenesis. |
| 組織特異性 | Expressed at high levels in the heart, pancreas, adrenal gland and ovary and at low levels in placenta, liver, kidney, prostate, testis, small intestine, spleen and colon. In the kidney, expressed by the visceral epithelial cells of the glomeruli. A widespread expression is also seen in the medial smooth muscle cells of arteries and arterioles, as well as in smooth muscle cells of vasa rectae in the medullary area. Expressed in the adventitial connective tissue surrounding the suprarenal artery. In chronic obstructive nephropathy, a persistent expression is seen in glomerular visceral epithelial cells and vascular smooth muscle cells, as well as de novo expression by periglomerular interstitial cells and by some neointimal cells of arteriosclerotic vessels. Expression in normal prostate is seen preferentially in the mesenchyme of the gland while expression is increased and more profuse in prostate carcinoma. Expressed in many ovarian, lung, renal and brain cancer-derived cell lines. |
| 配列類似性 | Belongs to the PDGF/VEGF growth factor family. Contains 1 CUB domain. |
| 発生段階 | Not detectable in the earliest stages of glomerulogenesis, and not detected in the metanephric blastema or surrounding cortical interstitial cells. In later stages of glomerulogenesis, localized to epithelial cells transitioning from the early developing nephrons of the comma- and S-shaped stages to the visceral epithelial cells of differentiated glomeruli. In the developing pelvis, expressed at the basement membrane of immature collecting ducts and by presumptive fibroblastic cells in the interstitium. |
| 翻訳後修飾 | Proteolytic removal of the N-terminal CUB domain releasing the core domain is necessary for unmasking the receptor-binding epitopes of the core domain. Cleavage after Arg-247 or Arg-249 by urokinase plasminogen activator gives rise to the active form. |
| 細胞内局在 | Secreted. |

画像



All lanes : Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] (ab181845) at 1/1000 dilution

Lane 1 : Human fetal kidney lysate

Lane 2 : Human fetal heart lysate


Lane 3 : Human small intestine lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab136636](#)) at 1/500 dilution (Goat Anti-Rabbit IgG HRP)


Predicted band size: 43 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] (ab181845)

Immunohistochemical analysis of paraffin-embedded Mouse heart tissue labeling SCDGFB/PDGF-D using ab181845 at 1/100 dilution. A ready to use HRP Polymer for Rabbit IgG prediluted was used as secondary. Counterstain: Hematoxylin.


Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] (ab181845)

Immunohistochemical analysis of paraffin-embedded Human squamous cell lung carcinoma tissue labeling SCDGFB/PDGF-D using ab181845 at 1/100 dilution. A ready to use HRP Polymer for Rabbit IgG prediluted was used as secondary. Counterstain: Hematoxylin.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.



Western blot - Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] (ab181845)


Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] (ab181845) at 1/1000 dilution + PC3 cell lysate at 10 µg

Secondary

Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 43 kDa

Why choose a recombinant antibody?



Anti-SCDGFB/PDGF-D antibody [EPR7266(2)] (ab181845)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors