


Anti-SDHB antibody [21A11AE7] ab14714

KO 評価済

★★★★★ [12 Abreviews](#) [252 References](#) [画像数 6](#)

製品の概要

製品名	Anti-SDHB antibody [21A11AE7]
製品の詳細	Mouse monoclonal [21A11AE7] to SDHB
由来種	Mouse
アプリケーション	適用あり: Flow Cyt, IHC-Fr, WB, ICC/IF
種交差性	交差種: Mouse, Rat, Cow, Human 交差が予測される動物種: Hamster, Pig, Zebrafish 
免疫原	Full length protein corresponding to Cow SDHB. Purified SDH from bovine heart.
ポジティブ・コントロール	WB: Human, bovine, rat, and mouse heart mitochondria; Isolated mitochondria from HepG2 cells. IHC-Fr: Normal ageing human colon tissue. Flow Cyt: HEK-293 cells. ICC/IF: HeLa cells. HEK293 cells (HEK293-SDHB KO used as a negative cell line)
特記事項	<p>This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact orders@abcam.com.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> <p>Product was previously marketed under the MitoSciences sub-brand.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.5 Preservative: 0.02% Sodium azide Constituent: HEPES buffered saline

精製度	IgG fraction
特記事項 (精製)	Near homogeneity as judged by SDS-PAGE. The antibody was produced in vitro using hybridomas grown in serum-free medium, and then purified by biochemical fractionation.
ポリ/モノ	モノクローナル
クローン名	21A11AE7
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab14714の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 1µg for 10 ⁶ cells. ab170191 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.
IHC-Fr		Use at an assay dependent concentration.
WB	★★★★★ (9)	Use a concentration of 5 µg/ml. Detects a band of approximately 28 kDa.
ICC/IF	★★★★★ (1)	Use a concentration of 1 µg/ml.

ターゲット情報

機能	Iron-sulfur protein (IP) subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone (coenzyme Q).
パスウェイ	Carbohydrate metabolism; tricarboxylic acid cycle; fumarate from succinate (eukaryal route): step 1/1.
関連疾患	<p>Defects in SDHB are a cause of susceptibility to pheochromocytoma (PCC) [MIM:171300]. A catecholamine-producing tumor of chromaffin tissue of the adrenal medulla or sympathetic paraganglia. The cardinal symptom, reflecting the increased secretion of epinephrine and norepinephrine, is hypertension, which may be persistent or intermittent.</p> <p>Defects in SDHB are the cause of hereditary paragangliomas type 4 (PGL4) [MIM:115310]; also known as familial non-chromaffin paragangliomas type 4. Paragangliomas refer to rare and mostly benign tumors that arise from any component of the neuroendocrine system. PGL4 is characterized by the development of mostly benign, highly vascular, slow growing tumors in the head and neck. In the head and neck region, the carotid body is the largest of all paraganglia and is also the most common site of the tumors.</p> <p>Defects in SDHB are a cause of paraganglioma and gastric stromal sarcoma (PGGSS) [MIM:606864]; also called Carney-Stratakis syndrome. Gastrointestinal stromal tumors may be sporadic or inherited in an autosomal dominant manner, alone or as a component of a syndrome associated with other tumors, such as in the context of neurofibromatosis type 1 (NF1). Patients</p>

have both gastrointestinal stromal tumors and paragangliomas. Susceptibility to the tumors was inherited in an apparently autosomal dominant manner, with incomplete penetrance. Defects in SDHB are a cause of Cowden-like syndrome (CWDLS) [MIM:612359]. Cowden-like syndrome is a cancer predisposition syndrome associated with elevated risk for tumors of the breast, thyroid, kidney and uterus.


配列類似性

Belongs to the succinate dehydrogenase/fumarate reductase iron-sulfur protein family. Contains 1 2Fe-2S ferredoxin-type domain. Contains 1 4Fe-4S ferredoxin-type domain.

細胞内局在


Mitochondrion inner membrane.

画像



ab14714 staining SDHB in wild-type HEK293 cells (top panel) and SDHB knockout HEK293 cells (bottom panel). The cells were fixed with 100% methanol (5 min), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated with ab14714 at 5µg concentration and [ab6046](#) (Rabbit polyclonal to beta Tubulin) at 1/1000 dilution overnight at +4°C, followed by a further incubation at room temperature for 1h with a goat secondary antibody to mouse IgG (Alexa Fluor® 488) ([ab150117](#)) at 2 µg/ml (shown in green) and a goat secondary antibody to rabbit IgG (Alexa Fluor® 594) ([ab150080](#)) at 2 µg/ml (shown in red). Nuclear DNA was labelled in blue with DAPI.

Image was taken with a confocal microscope (Leica-Microsystems, TCS SP8).



All lanes : Anti-SDHB antibody [21A11AE7] (ab14714) at 5 µg/ml

Lane 1 : Isolated mitochondria from Human heart at 5 µg

Lane 2 : Isolated mitochondria from Bovine Heart at 1 µg

Lane 3 : Isolated mitochondria from Rat heart at 10 µg

Lane 4 : Isolated mitochondria from Mouse heart at 10 µg

Lane 5 : Isolated mitochondria from HepG2 (human liver hepatocellular carcinoma cell line) cells at 20 µg


Secondary

All lanes : Goat anti-Mouse secondary

Observed band size: 28 kDa


Additional bands at: 55 kDa. We are unsure as to the identity of these extra bands.

Extra bands in the mouse sample (lane 4) are due to the reaction of the IgG-specific goat anti-mouse secondary antibody with residual mouse blood in the heart tissue, as it is very difficult to entirely remove the blood from these small organs.




Immunohistochemistry (Frozen sections) - Anti-SDHB antibody [21A11AE7] (ab14714)

ab14714 staining SDHB in normal ageing human colon tissue by Immunohistochemistry (Frozen sections).




Immunocytochemistry/ Immunofluorescence - Anti-SDHB antibody [21A11AE7] (ab14714)

Mitochondrial localization of complex II visualized by immunocytochemistry using anti-complex II subunit 30 kDa Ip mAb 21A11 (ab14714). Cells were fixed, permeabilized and then labeled with ab14714 followed by an Alexa Fluor® 488-conjugated-goat-anti-mouse IgG2a isotype specific secondary antibody.



Flow Cytometry - Anti-SDHB antibody [21A11AE7] (ab14714)

Overlay histogram showing HEK-293 (human epithelial cell line from embryonic kidney) cells stained with ab14714 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab14714, 1µg/1x10⁶ cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) (**ab96879**) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG2a [ICIGG2A] (**ab91361**, 1µg/1x10⁶ cells) used under the same conditions. Acquisition of >5,000 events was performed.



Immunocytochemistry/ Immunofluorescence - Anti-SDHB antibody [21A11AE7] (ab14714)

ICC/IF image of ab14714 stained HeLa (human epithelial cell line from cervix adenocarcinoma) cells. The cells were 4% formaldehyde fixed (10 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (ab14714, 1µg/ml) overnight at +4°C. The secondary antibody (green) was Alexa Fluor® 488 goat anti-mouse IgG (H+L) used at a 1/1000 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors