

Anti-SFXN1/TCC antibody ab127751

画像検索

医薬用外劇物

製品の概要

製品名	Anti-SFXN1/TCC antibody
製品の詳細	Rabbit polyclonal to SFXN1/TCC
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment corresponding to Human SFXN1/TCC aa 1-146. Database link: NP_073591
ポジティブ・コントロール	293T whole cell lysate, A431 whole cell lysate, Jurkat whole cell lysate, Raji whole cell lysate, Human hepatoma tissue.
特記事項	Keep as concentrated solution. The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 20% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab127751の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 36 kDa.
IHC-P		1/100 - 1/500.

ターゲット情報

機能	Might be involved in the transport of a component required for iron utilization into or out of the mitochondria.
配列類似性	Belongs to the sideroflexin family.
細胞内局在	Mitochondrion membrane.


画像



Anti-SFXN1/TCC antibody (ab127751) at 1/1000 dilution + 293T whole cell lysate at 30 µg

Predicted band size: 36 kDa

10% SDS PAGE



ab127751, at 1/500, staining SFXN1/TCC in paraffin embedded Human hepatoma tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SFXN1/TCC antibody (ab127751)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors