


Anti-SH3BGRL3 antibody ab122138

画像数 3

製品の概要

製品名	Anti-SH3BGRL3 antibody
製品の詳細	Rabbit polyclonal to SH3BGRL3
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB, ICC/IF
種交差性	交差種: Human
免疫原	antigen sequence: SEVTRILDGK RIQYQLVDIS QDNALRDEMRL ALAGNPKATP PQIVNGDQYC GDYELFVEAV EQNTLQEFLK LA, corresponding to amino acids 22-93 of Human SH3BGRL3  Run BLAST with   Run BLAST with 
ポジティブ・コントロール	Human smooth muscle tissue; RT4 cell lysates and Human liver and tonsil tissue lysates; U2OS cells
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab122138の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/600. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use a concentration of 0.04 - 0.4 µg/ml.
ICC/IF		Use a concentration of 1 - 4 µg/ml. Fixed/permeabilized with PFA/Triton X-100

ターゲット情報

機能	Could act as a modulator of glutaredoxin biological activity.
組織特異性	Ubiquitous.
配列類似性	Belongs to the SH3BGR family. Contains 1 glutaredoxin domain.
細胞内局在	Cytoplasm. Nucleus.

画像



All lanes : Anti-SH3BGR3 antibody (ab122138) at 1/250 dilution


Lane 1 : RT4 cell lysates

Lane 2 : U251MG cell lysates

Lane 3 : Human plasma


Lane 4 : Human liver tissue lysate

Lane 5 : Human tonsil tissue lysate



ab122138, at 1/600 dilution, staining SH3BGRL3 in Human smooth muscle tissue by Immunohistochemistry using paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SH3BGRL3 antibody (ab122138)



ab122138 at 4 µg/ml staining SH3BGRL3 in PFA/Triton X-100 fixed/permeabilized U2OS cells by Immunofluorescence. Antibody staining is shown in green

Immunocytochemistry/ Immunofluorescence - Anti-SH3BGRL3 antibody (ab122138)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors