


Anti-SHP1 (phospho S591) antibody ab41436

1 References [画像数 1](#)

製品の概要

製品名	Anti-SHP1 (phospho S591) antibody
製品の詳細	Rabbit polyclonal to SHP1 (phospho S591)
由来種	Rabbit
特異性	Phospho-Ser591-SHP1 specific. This antibody detects a 68 kDa band on SDS-PAGE immunoblots of human Jurkat and A431 cells treated with pervanadate, as well as A431 cells treated with EGF.
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Synthetic peptide corresponding to Human SHP1 (phospho S591). Database link: P29350
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium azide Constituents: PBS, 50% Glycerol, 0.1% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was cross-adsorbed to a non specific phospho-serine peptide then affinity purified using phospho-Ser591-SHP1 peptide (immunogen).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab41436の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
WB		1/1000. Detects a band of approximately 68 kDa (predicted molecular weight: 68 kDa).

ターゲット情報

機能	Plays a key role in hematopoiesis. This PTPase activity may directly link growth factor receptors and other signaling proteins through protein-tyrosine phosphorylation. The SH2 regions may interact with other cellular components to modulate its own phosphatase activity against interacting substrates. Together with MTUS1, induces UBE2V2 expression upon angiotensin II stimulation.
組織特異性	Isoform 1 is expressed in hematopoietic cells. Isoform 2 is expressed in non-hematopoietic cells.
配列類似性	Belongs to the protein-tyrosine phosphatase family. Non-receptor class 2 subfamily. Contains 2 SH2 domains. Contains 1 tyrosine-protein phosphatase domain.
翻訳後修飾	Phosphorylated on serine and tyrosine residues.
細胞内局在	Cytoplasm. Nucleus. In neurons, translocates into the nucleus after treatment with angiotensin II.

画像



Lanes 1-2 : anti-SHP1 (C-terminal)

Lanes 3-6 : Anti-SHP1 (phospho S591) antibody (ab41436) at 1/1000 dilution

Lanes 1 & 3 : human Jurkat cells treated with 1 mM pervanadate for 30 min.

Lanes 2 & 4 : human Jurkat cells treated with 1 mM pervanadate for 30 min, blot exposed to lambda phosphatase before probing

Lane 5 : human Jurkat cells treated with 1 mM pervanadate for 30 min. with phospho-Ser591-SHP1 peptide

Lane 6 : human Jurkat cells treated with 1 mM pervanadate for 30 min. with non specific phospho-serine peptide

Predicted band size: 68 kDa

Observed band size: 68 kDa

Membrane was incubated with diluted antibody in 5% non-fat milk, PBS, 0.04% Tween20 for 1 hour at room temperature.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors