

# Anti-Streptolysin antibody [6D11] ab23501

## 4 References

### 製品の概要

製品名	Anti-Streptolysin antibody [6D11]
製品の詳細	Mouse monoclonal [6D11] to Streptolysin
由来種	Mouse
特異性	ab23501 is specific for oxidized streptolysin but cross reacts with untreated streptolysin
アプリケーション	<b>適用あり:</b> ELISA
種交差性	<b>交差種:</b> Streptococcus
免疫原	Cysteine chloride oxidized streptolysin.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: 0.0268% PBS, 2.9% Sodium chloride</p>
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	6D11
ミエローマ	x63-Ag8.653
アイソタイプ	IgG1
軽鎖の種類	kappa

## アプリケーション

**The Abpromise guarantee**      Abpromise保証は、次のテスト済みアプリケーションにおけるab23501の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/400.

## ターゲット情報

**関連性**      Bacteria of the Streptococcus genera produce the exotoxin streptolysin. Anti-streptolysin antibody titers in patients are measured as a marker of inflammation. Streptolysin is also an important tool in cell biology due to its hemolytic activity and the ability to permeabilize plasma membranes. Sulfhydryl-activated toxin. Is able to lyse cholesterol containing membranes. Can be reversibly inactivated by oxidation. Cholesterol is the receptor for the binding of these toxins to eukaryotic cell membranes.

**細胞内局在**      Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors