

### Anti-Syndecan 3 antibody ab63932

★ ★ ★ ★ ★ 1 Abreviews 7 References 画像数 1

医薬用外劇物

#### 製品の概要

製品名	Anti-Syndecan 3 antibody
製品の詳細	Rabbit polyclonal to Syndecan 3
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Rat
免疫原	Synthetic peptide surrounding amino acid 427 of human Syndecan 3
特記事項	<p>This product is manufactured by BioVision, an Abcam company and was previously called 3824 Syndecan-3 Antibody. 3824-100 is the same size as the 100 µg size of ab63932.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.01% Thimerosal (merthiolate)</p> <p>Constituents: PBS, 30% Glycerol (glycerin, glycerine), 0.5% BSA</p>
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

## The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab63932の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 - 4 µg/ml. Detects a band of approximately 46 kDa (predicted molecular weight: 46 kDa).

## ターゲット情報

機能	Cell surface proteoglycan that may bear heparan sulfate (By similarity). May have a role in the organization of cell shape by affecting the actin cytoskeleton, possibly by transferring signals from the cell surface in a sugar-dependent mechanism.
組織特異性	Expressed in the nervous system, the adrenal gland, and the spleen.
配列類似性	Belongs to the syndecan proteoglycan family.
翻訳後修飾	O-glycosylated within the Thr/Ser-rich region which could interact with lectin domains on other molecules.
細胞内局在	Membrane.

## 画像



Anti-Syndecan 3 antibody (ab63932) at 1 µg/ml + Rat kidney tissue lysate at 30 µg

**Predicted band size:** 46 kDa

**Observed band size:** 46 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors