

Anti-Syntaxin 8 antibody ab102798

画像数 2

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-Syntaxin 8 antibody |
| 製品の詳細 | Rabbit polyclonal to Syntaxin 8 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Human |
| 免疫原 | Recombinant full length protein within Human Syntaxin 8 aa 1-250. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please contact our Scientific Support team to discuss your requirements. Run BLAST with Run BLAST with |
| ポジティブ・コントロール | Human kidney tissue lysate. Transfected 293T cell line. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | pH: 7.40 Constituent: 100% PBS |
| 精製度 | Protein A purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab102798の使用に適用されず
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 27 kDa. |

ターゲット情報

| | |
|-------|---|
| 機能 | Vesicle trafficking protein that functions in the early secretory pathway, possibly by mediating retrograde transport from cis-Golgi membranes to the ER. |
| 組織特異性 | Highly expressed in heart. Also found in brain, kidney, liver, lung, placenta, skeletal muscle, spleen and pancreas. |
| 配列類似性 | Belongs to the syntaxin family. Contains 1 t-SNARE coiled-coil homology domain. |
| 翻訳後修飾 | Ubiquitinated by HECTD3. |
| 細胞内局在 | Membrane. Preferentially associated with the early endosome. To lesser extends, also present in late endosome, the plasma membrane and coated pits. |


画像



Western blot - Anti-Syntaxin 8 antibody (ab102798)

Anti-Syntaxin 8 antibody (ab102798) at 5 µg/ml + Human kidney tissue lysate at 50 µg

Predicted band size: 27 kDa



Western blot - Anti-Syntaxin 8 antibody (ab102798)

All lanes : Anti-Syntaxin 8 antibody (ab102798) at 5 µg/ml

Lane 1 : Syntaxin 8 transfected 293T lysate

Lane 2 : Non-transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 27 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors