

### Anti-TAF7 antibody [2C5] ab57494

★★★★★ [2 Abreviews](#) [3 References](#) [画像数 5](#)

#### 製品の概要

製品名	Anti-TAF7 antibody [2C5]
製品の詳細	Mouse monoclonal [2C5] to TAF7
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB, IHC-P, Flow Cyt
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment: FIWNHGITLP LKNVRKRRFR KTAKKKYIES PDVEKEVKRL LSTDAEAVST RWEIIAEDT KEAENQGLDI SSPGMSGHRQ GHDSLEHDEL REIFN, corresponding to amino acids 130-224 of Human TAF7 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

#### 特記事項

This product was changed from ascites to tissue culture supernatant on 12<sup>th</sup> June 2019. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.4
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	2C5
アイソタイプ	IgG2a

アプリケーション


The Abpromise guarantee      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab57494の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use at an assay dependent concentration. Predicted molecular weight: 40 kDa.
IHC-P	★★★★★ (1)	Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. <b>ab170191</b> - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.


ターゲット情報

機能	Functions as a component of the DNA-binding general transcription factor complex TFIID, a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. Present in both of the previously described TFIID species which either lack or contain TAFII30 (TFIID alpha and TFIID beta respectively).
組織特異性	Ubiquitous.
配列類似性	Belongs to the TAF7 family.
ドメイン	The [KR]-[STA]-K motif is specifically recognized by the SETD7 methyltransferase, which methylates Lys-5 in vitro.
細胞内局在	Nucleus.

画像



TAF7 antibody (ab57494) at 1ug/lane + MCF-7 cell lysate at 25ug/lane.  
This image was generated using the ascites version of the product.




Flow Cytometry - Anti-TAF7 antibody [2C5]  
(ab57494)

Overlay histogram showing HeLa cells stained with ab57494 (red line). The cells were fixed with 4% paraformaldehyde (10 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab57494, 1µg/1x10<sup>6</sup> cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG2a [ICIGG2A] ([ab91361](#), 1µg/1x10<sup>6</sup> cells) used under the same conditions.

Unlabelled sample (blue line) was also used as a control. Acquisition of >5,000 events were collected using a 20mW Argon ion laser (488nm) and 525/30 bandpass filter.


This image was generated using the ascites version of the product.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TAF7 antibody [2C5]  
(ab57494)

TAF7 antibody (ab57494) used in immunohistochemistry at 1ug/ml on formalin fixed and paraffin embedded human lymph node.

This image was generated using the ascites version of the product.



Western blot - Anti-TAF7 antibody [2C5] (ab57494)

Image kindly supplied by Dr Karin Birkenkamp-Demtroeder through Abreview

**Lane 2 :** Anti-TAF7 antibody [2C5] (ab57494) at 1/1250 dilution

**Lane 3 :** Anti-TAF7 antibody [2C5] (ab57494) at 1/2500 dilution

**Lane 4 :** Anti-TAF7 antibody [2C5] (ab57494) at 1/5000 dilution

**Lane 5 :** Anti-TAF7 antibody [2C5] (ab57494) at 1/10000 dilution

**All lanes :** whole cell lysate prepared from SW780 badder cancer cells

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat anti-mouse IgG conjugated to HRP

Developed using the ECL technique.

Performed under reducing conditions.

**Predicted band size:** 40 kDa


**Observed band size:** 40 kDa

**Additional bands at:** 90 kDa. We are unsure as to the identity of these extra bands.

**Exposure time:** 10 minutes

Lane 1: All Blue Biorad. Gel run under denaturing conditions with a 4-12% gradient. Blocking step performed using 5% milk for one hour at room temperature. Primary antibody incubated for 16 hours.

This image was generated using the ascites version of the product.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TAF7 antibody [2C5] (ab57494)

This image was kindly supplied by Dr Karin Birkenkamp-Demtroeder by Abreview

ab57494 staining TAF7 in human bladder cancer tissue sections by Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections). Tissue was fixed with formaldehyde and a heat mediated antigen retrieval step was performed using TEG. Samples were then blocked with 1% serum for 1 hour at room temperature followed by incubation with the primary antibody at a 1/3500 dilution for 1 hour. An undiluted HRP-conjugated goat anti-mouse IgG polyclonal was used as secondary antibody.

This image was generated using the ascites version of the product.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors