

Anti-Thrombospondin 2 antibody [MM0568-9T19] ab89805

[1 Abreviews](#) [1 References](#) [画像数 1](#)

製品の概要

製品名	Anti-Thrombospondin 2 antibody [MM0568-9T19]
製品の詳細	Mouse monoclonal [MM0568-9T19] to Thrombospondin 2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant full length Human Thrombospondin 2
ポジティブ・コントロール	WB: Human placenta tissue lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
精製度	Protein G purified
特記事項 (精製)	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and filtered through a 0.2 µm filter.
ポリ/モノ	モノクローナル
クローン名	MM0568-9T19
アイソタイプ	IgG2

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab89805の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 130 kDa.

ターゲット情報

機能	Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. Can bind to fibrinogen, fibronectin, laminin and type V collagen.
組織特異性	High expression in intervertebral disk tissue.
関連疾患	Genetic variations in THBS2 may be a cause of susceptibility to intervertebral disc disease (IDD) [MIM:603932]; also known as lumbar disk herniation (LDH). IDD is one of the most common musculo-skeletal disorders and the predominant cause of low-back pain and unilateral leg pain.
配列類似性	Belongs to the thrombospondin family. Contains 3 EGF-like domains. Contains 1 TSP C-terminal (TSPC) domain. Contains 1 TSP N-terminal (TSPN) domain. Contains 3 TSP type-1 domains. Contains 8 TSP type-3 repeats. Contains 1 VWFC domain.

画像



Anti-Thrombospondin 2 antibody [MM0568-9T19] (ab89805) +
Human placenta tissue lysate

Predicted band size: 130 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors