

Anti-TIGIT antibody [EPR25804-63] - BSA and Azide free (Capture) ab290682

リコンビナント RabMAb

画像数 2

製品の概要

| | |
|----------|--|
| 製品名 | Anti-TIGIT antibody [EPR25804-63] - BSA and Azide free (Capture) |
| 製品の詳細 | Rabbit monoclonal [EPR25804-63] to TIGIT - BSA and Azide free (Capture) |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: Sandwich ELISA |
| 種交差性 | 交差種: Human |
| 免疫原 | Recombinant fragment. This information is proprietary to Abcam and/or its suppliers. |
| 特記事項 | <p>ab290682 is a BSA and Azide Free antibody supplied in an unconjugated format and it is suitable for sandwich ELISAs to quantify Human TIGIT. The recommended pair for sandwich ELISA is:</p> <p>Capture: ab290682 , Human TIGIT Capture Antibody (unconjugated) Detector: ab2906823, Human TIGIT Detector Antibody (unconjugated)</p> <p>The reference range value is 5.86-375 pg/mL.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> |

製品の特性

| | |
|----------|----------------------------------|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | pH: 7.2 Constituent: 100% PBS |
| キャリア・フリー | はい |

| | |
|--------|--------------------|
| 精製度 | Protein A purified |
| ポリモノ | モノクローナル |
| クローン名 | EPR25804-63 |
| アイソタイプ | IgG |

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab290682の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------------|-----------|--|
| Sandwich ELISA | | Use at an assay dependent concentration. Can be paired for Sandwich ELISA with PE / R-Phycoerythrin Conjugation Kit - Lightning-Link® (ab102918) and Streptavidin Conjugation Kit - Lightning-Link® (ab102921) and Biotinylation Kit / Biotin Conjugation Kit (Fast and Type A) - Lightning-Link® (ab201795) and Rabbit monoclonal [EPR25804-72] to TIGIT - BSA and Azide free (Detector) (ab290683) . |

ターゲット情報

| | |
|-------|---|
| 機能 | Binds with high affinity to the poliovirus receptor (PVR) which causes increased secretion of IL10 and decreased secretion of IL12B and suppresses T cell activation by promoting the generation of mature immunoregulatory dendritic cells. |
| 組織特異性 | Expressed at low levels on peripheral memory and regulatory CD4+ T cells and NK cells and is up-regulated following activation of these cells (at protein level). |
| 配列類似性 | Contains 1 Ig-like V-type (immunoglobulin-like) domain. |
| ドメイン | Contains 1 copy of a cytoplasmic motif that is referred to as the immunoreceptor tyrosine-based inhibitor motif (ITIM). This motif is involved in modulation of cellular responses. The phosphorylated ITIM motif can bind the SH2 domain of several SH2-containing phosphatases. |
| 細胞内局在 | Cell membrane. |

画像



Sandwich ELISA of **ab290706** with the capture antibody dilution at 2 µg/mL and detector antibody dilution at 0.5 µg/mL.

Sandwich ELISA - Anti-TIGIT antibody [EPR25804-63] - BSA and Azide free (Capture) (ab290682)

Why choose a recombinant antibody?

| | |
|--|--|
|  <p>Research with confidence Consistent and reproducible results</p> |  <p>Long-term and scalable supply Recombinant technology</p> |
|  <p>Success from the first experiment Confirmed specificity</p> |  <p>Ethical standards compliant Animal-free production</p> |

Anti-TIGIT antibody [EPR25804-63] - BSA and Azide free (Capture) (ab290682)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors