

# Anti-TMEM143 antibody ab121431

画像数 3

### 製品の概要

製品名	Anti-TMEM143 antibody
製品の詳細	Rabbit polyclonal to TMEM143
由来種	Rabbit
アプリケーション	<b>適用あり:</b> ICC/IF, IHC-P, WB
種交差性	<b>交差種:</b> Human
免疫原	fragment, corresponding to internal sequence amino acids 57-197 of Human TMEM143, isoform 1.
ポジティブ・コントロール	Human colon tissue; Transfected HEK293T cells.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)</p>
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

The Abpromise guarantee

**Abpromise保証は、** 次のテスト済みアプリケーションにおけるab121431の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use a concentration of 0.04 - 0.4 µg/ml.

## ターゲット情報

細胞内局在


Membrane.

## 画像




Immunofluorescent staining of human cell line U-251 MG shows positivity in nucleoli.

Immunocytochemistry/ Immunofluorescence - Anti-TMEM143 antibody (ab121431)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMEM43 antibody (ab121431)

ab121431 at a 1/250 dilution staining TMEM43 in paraffin-embedded Human colon tissue, by Immunohistochemistry.



Western blot - Anti-TMEM43 antibody (ab121431)

**All lanes :** Anti-TMEM43 antibody (ab121431) at 1/250 dilution

**Lane 1 :** Vector only transfected HEK293T lysate

**Lane 2 :** TMEM43-transfected HEK293T over-expression Lysate  
(Co-expressed with a C-terminal myc-DDK tag ~3.1 kDa)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors