

Anti-TMEM223 antibody ab121434

画像数 2

製品の概要

製品名	Anti-TMEM223 antibody
製品の詳細	Rabbit polyclonal to TMEM223
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, IHC-P
種交差性	交差種: Human
免疫原	fragment corresponding to amino acids 119-202 (RSVRSVVLRA GGQQVLTTH APFGLGAHFT VPLKQVSCMA HRGEVPAMLP LKVKGRRFYF LLDKTGHFPN TKLFDNTVGA YRSL) of Human TMEM223.

[Run BLAST with Expasy](#) [Run BLAST with NCBI](#)

ポジティブ・コントロール IHC: Human duodenum tissue. ICC/IF: human cell line MCF7

特記事項 The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab121434の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 0.25 - 0.2 µg/ml.
IHC-P		Use at an assay dependent concentration.

ターゲット情報


配列類似性

Belongs to the TMEM223 family.

細胞内局在


Membrane.

画像



Immunocytochemistry/immunofluorescence analysis of human cell line MCF7 labelling TMEM223 with ab121434 at 0.25 µg/mL.

Immunocytochemistry/ Immunofluorescence - Anti-TMEM223 antibody (ab121434)



ab121434 staining TMEM223 in paraffin-embedded Human duodenum tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMEM223 antibody (ab121434)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors