

# Anti-TMX3 antibody ab121414

画像数 2

### 製品の概要

製品名	Anti-TMX3 antibody
製品の詳細	Rabbit polyclonal to TMX3
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, WB
種交差性	<b>交差種:</b> Human
免疫原	antigen sequence, corresponding to N terminal amino acids 16-143 of Human TMX3.
ポジティブ・コントロール	WB: RT 4, U 251 MG and Human tonsil lysates IHC: Human stomach tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab121414の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/10 - 1/20. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		Use a concentration of 0.04 - 0.4 µg/ml.

## ターゲット情報


機能	Probable disulfide isomerase, which participates in the folding of proteins containing disulfide bonds. May act as a dithiol oxidase.
組織特異性	Widely expressed. Expressed in brain, testis, lung, skin, kidney, uterus, bone, stomach, liver, prostate, placenta, eye and muscle.
配列類似性	Belongs to the protein disulfide isomerase family. Contains 1 thioredoxin domain.
ドメイン	The di-lysine motif confers endoplasmic reticulum localization for type I membrane proteins.
翻訳後修飾	N-glycosylated.
細胞内局在	Endoplasmic reticulum membrane.

## 画像



ab121414 at 1/10 dilution staining TMX3 in paraffin-embedded Human stomach tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMX3 antibody (ab121414)



**All lanes :** Anti-TMX3 antibody (ab121414) at 1/250 dilution

**Lane 1 :** RT 4 cell lysate

**Lane 2 :** U 251 MG cell lysate

**Lane 3 :** Human Plasma lysate

**Lane 4 :** Human liver tissue lysate

**Lane 5 :** Human tonsil tissue lysate

Developed using the ECL technique.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors