

Anti-USP37 antibody ab190184

[4 References](#) [画像数 2](#)

製品の概要

製品名	Anti-USP37 antibody
製品の詳細	Rabbit polyclonal to USP37
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat, Human, Xenopus laevis
免疫原	Recombinant fragment corresponding to Human USP37 aa 1-300. Database link: Q86T82 Run BLAST with Run BLAST with
ポジティブ・コントロール	U2OS cells transfected with USP37 siRNA lysate, Mouse embryo fibroblasts cell lysate and INS1 cell lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium azide
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab190184の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/2000. Predicted molecular weight: 110 kDa.

ターゲット情報

組織特異性


Expressed in brain and prostate.

配列類似性

Belongs to the peptidase C19 family.

Contains 3 UIM (ubiquitin-interacting motif) repeats.

画像



Western blot - Anti-USP37 antibody (ab190184)

All lanes : Anti-USP37 antibody (ab190184) at 1/2000 dilution

Lane 1 : Human U2OS cell lysate si Control


Lane 2 : Human U2OS cells transfected with USP37 siRNA

Lane 3 : Mouse Embryo Fibroblaste cell Lysate

Lane 4 : Rat INS1 cell lysate

Lysates/proteins at 40 µg per lane.

Predicted band size: 110 kDa



Western blot - Anti-USP37 antibody (ab190184)

All lanes : Anti-USP37 antibody (ab190184) at 1/2000 dilution

Lane 1 : Xenopus laevis crude egg extracts

Lysates/proteins at 50 µg per lane.

Predicted band size: 110 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors