




Anti-V5 tag antibody [SV5-P-K] ab206557

リコンビナント

画像数 2

製品の概要

製品名	Anti-V5 tag antibody [SV5-P-K]
製品の詳細	Mouse monoclonal [SV5-P-K] to V5 tag
由来種	Mouse
アプリケーション	適用あり: ELISA
種交差性	交差種: Species independent
免疫原	Synthetic peptide corresponding to V5 tag. (Paramyxovirus simian virus 5).  Run BLAST with  Run BLAST with 
特記事項	ab206557 is a recombinant monoclonal antibody to V5 epitope tag, manufactured using variable regions (i.e. specificity) from the hybridoma SV5-P-K.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Proclin 300 Constituent: 99% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	SV5-P-K
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab206557の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.

ターゲット情報

関連性

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPNPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

画像



ELISA - Anti-V5 tag antibody [SV5-P-K] (ab206557)

ab206557 loaded at 1000 ng/ml in first well and an 11 point serial dilution performed in duplicate for each antibody with last well left blank. Signal detection was performed using HRP-conjugated goat anti-mouse. Plate coated with 5ug/ml of Multitag protein.

Why choose a recombinant antibody?

Research with confidence
Consistent and reproducible results

Long-term and scalable supply
Recombinant technology

Success from the first experiment
Confirmed specificity

Ethical standards compliant
Animal-free production

Anti-V5 tag antibody [SV5-P-K] (ab206557)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors