

# Anti-YB1 antibody [EPR22677-250] - BSA and Azide free (Capture) ab270328

リコンビナント RabMAb

画像数 2

### 製品の概要

製品名	Anti-YB1 antibody [EPR22677-250] - BSA and Azide free (Capture)
製品の詳細	Rabbit monoclonal [EPR22677-250] to YB1 - BSA and Azide free (Capture)
由来種	Rabbit
アプリケーション	<b>適用あり:</b> Sandwich ELISA
種交差性	<b>交差種:</b> Human
免疫原	Recombinant full length protein. This information is proprietary to Abcam and/or its suppliers.
特記事項	<p>ab270328 is a BSA and Azide Free antibody supplied in an unconjugated format and it is suitable for sandwich ELISAs to quantify Human YB1. The recommended pair for sandwich ELISA is:</p> <p>Capture: ab270328, Human YB1 Capture Antibody (unconjugated) Detector: <b>ab270329</b>, Human YB1 Detector Antibody (unconjugated)</p> <p>The reference range value is 0.9375 - 60 pg/ml.</p> <p>Our <b>carrier-free</b> antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our <b>conjugation kits</b> for antibody conjugates that are ready-to-use in as little as 20 minutes with &lt;1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>The recommended antibody orientation is based on internal optimization for ELISA-based assays. Antibody orientation is assay dependent and needs to be optimized for each assay type. Please note that the range provided for this antibody is only an estimation based on the performance of the product using the recommended antibody pair. Performance of the antibody pair will depend on the specific characteristics of your assay. We guarantee the product works in sandwich ELISA, but we do not guarantee the sensitivity or dynamic range of the antibody in your assay.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p>

- High batch-to-batch consistency and reproducibility
- Improved sensitivity and specificity
- Long-term security of supply
- Animal-free production

For more information [see here](#).

Our RabMAb<sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb<sup>®</sup> patents](#).

## 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Constituent: 100% PBS
キャリア・フリー	はい
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR22677-250
アイソタイプ	IgG

## アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab270328の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use at an assay dependent concentration. Can be paired for Sandwich ELISA with <a href="#">PE / R-Phycoerythrin Conjugation Kit - Lightning-Link<sup>®</sup> (ab102918)</a> and <a href="#">Biotinylation Kit / Biotin Conjugation Kit (Fast and Type A) - Lightning-Link<sup>®</sup> (ab201795)</a> .

## ターゲット情報

**機能**

Mediates pre-mRNA alternative splicing regulation. Binds to splice sites in pre-mRNA and regulates splice site selection. Binds and stabilizes cytoplasmic mRNA. Contributes to the regulation of translation by modulating the interaction between the mRNA and eukaryotic initiation factors (By similarity). Regulates the transcription of numerous genes. Its transcriptional activity on the multidrug resistance gene MDR1 is enhanced in presence of the APEX1 acetylated form at 'Lys-6' and 'Lys-7'. Binds to promoters that contain a Y-box (5'-CTGATTGGCCAA-3'), such as MDR1 and HLA class II genes. Promotes separation of DNA strands that contain mismatches or are modified by cisplatin. Has endonucleolytic activity and can introduce nicks or breaks into double-stranded DNA (in vitro). May play a role in DNA repair. Component of the CRD-mediated complex that promotes MYC mRNA stability.

The secreted form acts as an extracellular mitogen and stimulates cell migration and proliferation.

## 配列類似性

Contains 1 CSD (cold-shock) domain.

## 翻訳後修飾

Ubiquitinated by RBBP6; leading to a decrease of YBX1 transactivational ability.


In the absence of phosphorylation the protein is retained in the cytoplasm.

Cleaved by a 20S proteasomal protease in response to agents that damage DNA. Cleavage takes place in the absence of ubiquitination and ATP. The resulting N-terminal fragment accumulates in the nucleus.

## 細胞内局在

Cytoplasm. Nucleus. Cytoplasmic granule. Secreted. Localized in cytoplasmic mRNP granules containing untranslated mRNAs. Shuttles between nucleus and cytoplasm. Predominantly cytoplasmic in proliferating cells. Cytotoxic stress and DNA damage enhance translocation to the nucleus. Localized with DDX1, MBNL1 and TIAL1 in stress granules upon stress. Secreted by mesangial and monocytic cells after inflammatory challenges. Translocates from the cytoplasm to the nucleus after and colocalizes with APEX1 in nuclear speckles after genotoxic stress.





## 画像



Representative standard curve from corresponding SimpleStep ELISA® Kit ([ab269544](#)).

Sandwich ELISA - Anti-YB1 antibody [EPR22677-250] - BSA and Azide free (Capture) (ab270328)

Why choose a recombinant antibody?

- **Research with confidence**  
Consistent and reproducible results
- **Long-term and scalable supply**  
Recombinant technology
- **Success from the first experiment**  
Confirmed specificity
- **Ethical standards compliant**  
Animal-free production

Anti-YB1 antibody [EPR22677-250] - BSA and Azide free (Capture) (ab270328)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors