

# Active Poly (Glu, Tyr) peptide ab204877

## 2 References

### 製品の詳細

製品名	Active Poly (Glu, Tyr) peptide
精製度	> 95 % Affinity purified. Purified by affinity chromatography
Animal free	No
由来	Synthetic
配列の追加情報	Molecular Weight: 20000 - 50000 Da

### 特性

Our **Abpromise guarantee** covers the use of **ab204877** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies
製品の状態	Liquid
備考	Ratio for molar concentration 4:1 Glu, Tyr (amount of each amino acid in the solution).

### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle. pH: 7.5 Constituent: 0.394% Tris HCl This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors