

D-VLK-p-NA peptide ab95075

製品の詳細

製品名	D-VLK-p-NA peptide
生理活性	Formula: H-D-Val-Leu-Lys-pNA. 2HCl
精製度	> 95 % SDS-PAGE.
Animal free	No
由来	Synthetic

特性

Our **Abpromise guarantee** covers the use of **ab95075** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	Formula: H-D-Val-Leu-Lys-pNA. 2HCl

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -80°C. Constituents: PBS, Tris buffered saline, 2.38% HEPES, 0.58% Sodium chloride
再構成	Resuspend to 10mM in 1mM HCl.

関連情報

関連性	D-VLK-p-NA is a chromogenic plasmin substrate
-----	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors