

IRS1 peptide ab204853

製品の詳細

製品名	IRS1 peptide
精製度	> 95 % HPLC.
Animal free	No
由来	Synthetic
配列	KKSRGDYMTMQIG
領域	979 to 989

特性

Our **Abpromise guarantee** covers the use of **ab204853** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies HPLC
製品の状態	Lyophilized
備考	ab204853 (IRS1 peptide) can be utilized as a substrate for the following active protein kinases:

- ab64298** (Active human Axl protein fragment)
- ab85755** (Active human Dystrophia myotonica protein kinase full length protein)
- ab69988** (Active Insulin Receptor R protein fragment)
- ab70687** (Active human Insulin Receptor protein fragment)
- ab70541** (Active human RON protein fragment)
- ab60337** (Active human Serine/threonine-protein kinase 4 full length protein)
- ab55710** (Active human LOK protein fragment)

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle.
再構成	Dilute peptide in distilled water to a final concentration of 1mg/ml

関連情報

機能	May mediate the control of various cellular processes by insulin. When phosphorylated by the insulin receptor binds specifically to various cellular proteins containing SH2 domains such as phosphatidylinositol 3-kinase p85 subunit or GRB2. Activates phosphatidylinositol 3-kinase when bound to the regulatory p85 subunit.
関連疾患	Polymorphisms in IRS1 may be involved in the etiology of non-insulin-dependent diabetes mellitus (NIDDM) [MIM:125853].
配列類似性	Contains 1 IRS-type PTB domain. Contains 1 PH domain.
翻訳後修飾	Serine phosphorylation of IRS1 is a mechanism for insulin resistance. Ser-312 phosphorylation inhibits insulin action through disruption of IRS1 interaction with the insulin receptor. Phosphorylation of Tyr-896 is required for GRB2-binding.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors