

Mouse CEBP Beta peptide ab39810

製品の詳細

製品名	Mouse CEBP Beta peptide
精製度	> 90 % SDS-PAGE.
Animal free	No
由来	Synthetic
生物種	Mouse
配列	MHRLLAWDAACLPPPAAFPKLH
領域	1 to 21

特性

Our **Abpromise guarantee** covers the use of **ab39810** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Blocking - Blocking peptide for Anti-CEBP Beta antibody [18F8] (ab15049)
	Neutralising
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle.
再構成	Reconstitute with 0.1 ml of distilled water.

関連情報

機能	Important transcriptional activator in the regulation of genes involved in immune and inflammatory responses. Specifically binds to an IL-1 response element in the IL-6 gene. NF-IL6 also binds to regulatory regions of several acute-phase and cytokines genes. It probably plays a role in the regulation of acute-phase reaction, inflammation and hemopoiesis. The consensus recognition site is 5'-T[TG]NNGNAA[TG]-3'. Functions in brown adipose tissue (BAT) differentiation.
組織特異性	Expressed at low levels in the lung, kidney and spleen.
配列類似性	Belongs to the bZIP family. C/EBP subfamily.

	Contains 1 bZIP domain.
ドメイン	the 9aaTAD motif is a transactivation domain present in a large number of yeast and animal transcription factors.
翻訳後修飾	Sumoylated by polymeric chains of SUMO2 or SUMO3.
細胞内局在	Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors