

Native Human Urokinase protein (FITC) ab93152

★★★★★ [1 Abreviews](#) [画像数 2](#)

製品の詳細

製品名	Native Human Urokinase protein (FITC)
精製度	> 95 % SDS-PAGE. > 98% pure by SDS PAGE
発現系	Native
タンパク質長	Full length protein
Animal free	No
由来	Native
生物種	Human
標識	FITC. Ex: 493nm, Em: 528nm

特性

Our **Abpromise guarantee** covers the use of **ab93152** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	Solubility: > 2 mg/mL Protect from light

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 5.00 Constituents: 0.41% Sodium acetate, 0.0292% EDTA, 0.58% Sodium chloride
------------	--

関連情報

機能	Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.
組織特異性	Expressed in the prostate gland and prostate cancers.
配列類似性	Belongs to the peptidase S1 family. Contains 1 EGF-like domain.

	Contains 1 kringle domain.
	Contains 1 peptidase S1 domain.
翻訳後修飾	Phosphorylation of Ser-158 and Ser-323 abolishes proadhesive ability but does not interfere with receptor binding.
細胞内局在	Secreted.


画像



10% SDS-PAGE

Lane 1: ab93152 (2 µg), non reduced

Lane 2: ab93152 (2 µg) plus Human PAI1 (stable mutant form) (10 µg), non reduced



10% SDS-PAGE

Lane 1: ab93152 (2 µg), non reduced

Lane 2: ab93152 (2 µg) plus Human PAI1 (stable mutant form) (10 µg), non reduced

Lane 3: Prestained standard

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors