

NLRC5 peptide ab196047

製品の詳細

製品名	NLRC5 peptide
Animal free	No
由来	Synthetic

特性

Our **Abpromise guarantee** covers the use of **ab196047** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Blocking - Blocking peptide for Anti-NLRC5 antibody (ab105411)
製品の状態	Liquid
備考	This is the blocking peptide for ab105411

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
------------	--

関連情報

機能	Probable regulator of the NF-kappa-B and type I interferon signaling pathways. May also regulate the type II interferon signaling pathway. Plays a role in homeostatic control of innate immunity and in antiviral defense mechanisms.
組織特異性	Expressed in spleen, thymus, lung, brain, tonsil, heart and prostate.
配列類似性	Belongs to the NLRP family. Contains 26 LRR (leucine-rich) repeats. Contains 1 NACHT domain.
細胞内局在	Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors