

Recombinant BK polyomavirus, strain AS, Major Capsid VP1 protein ab74567

2 References [画像数 1](#)

製品の詳細

製品名	Recombinant BK polyomavirus, strain AS, Major Capsid VP1 protein
精製度	> 90 % SDS-PAGE. Purified by ultracentrifugation.
発現系	Saccharomyces cerevisiae
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	BK virus
配列	<p>MAPTKRKGECPGAAPKKPKPEVQVPKLLIKGGVEVLEVKTGV DAITEVEC FLNPEMGDPDDNLRGYSQHLAENAFESDSPDRKMLPCYSTA RIPLPNLN EDLTCGNLLMWEAVTVKTEVIGITSMLNLHAGSQKVHENG KPVQGSNF HFFAVGGDPLEMQGVLMNYRTKYPQGTITPKNPTAQSVMNT DHKAYLDK NNAYPVECWIPDPSRNENTRYFGTYTGGENVPPVLHVTNTAT TVLLDEQG VGPLCKADSLYVSAADICGLFTNSSGTQQRGLARYFKIRLR KRSVKNPY PISFLLSDLINRRTQKVDGQPMYGMESQVEEVRVFDGTEQLP GDPDMIRY IDRQGQLQTKMV</p>

特性

Our **Abpromise guarantee** covers the use of **ab74567** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	Western blot
	SDS-PAGE

製品の状態

Lyophilized

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at 4°C prior to reconstitution. Store at -20°C long term. Reconstitute for long term storage.

Constituent: PBS

再構成

Reconstitute with deionized H₂O. With the addition of glycerol store at -20°C. After reconstitution can only be stored at 4°C for 2-3 days.

関連情報


関連性

BK polyomavirus was first isolated in 1971 in the urine of a renal transplant patient, initials B.K. BK polyomavirus, Major Capsid VP1 forms the icosahedral capsid of the virus protecting the viral genome. This capsid is composed of 360 VP1, assembled in 72 pentamers. Each pentamer are linked via the C-terminal end of VP1.

細胞内局在

Virion. Nucleus

画像



SDS-PAGE showing ab74567 at approximately 44kDa (3.2µg/lane).

Lane 1 shows the molecular weight ladder. From the bottom: 14.4, 18.4, 25.0, 35.0, 45.0, 66.2, 116.0 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors