

Recombinant Cytomegalovirus ICP36 protein ab43044

製品の詳細

製品名	Recombinant Cytomegalovirus ICP36 protein
生理活性	Immunoreactive with sera of CMV infected individuals, excellent antigen for detection of CMV with minimal specificity problems.
精製度	> 95 % SDS-PAGE. Sephadex derived purification.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
配列の追加情報	The protein contains the CMV ICP36 immunodominant regions.

特性

Our **Abpromise guarantee** covers the use of **ab43044** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	Western blot
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
	Constituents: 0.307% Glutathione, 0.79% Tris HCl, 50% Glycerol, 0.348% Sodium chloride

関連情報

関連性	In extracts from CMV-infected cells, the CMV DNA polymerase is found strongly associated with the additional polypeptide, ICP36. This protein has been identified as the CMV homolog of the herpes simplex virus type 1 UL42 gene product and may have a similar function. HCMV DNA polymerase and ICP36 expressed in the same system interact to form a stable complex. Moreover, ICP36 functions to stimulate the DNA polymerase activity in a template-dependent
-----	---

manner.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors