

Recombinant E. coli Ecotin protein ab92769

製品の詳細

製品名	Recombinant E. coli Ecotin protein
生理活性	Kinetic Data: 1:1 inhibition of trypsin (mole per mole)
精製度	> 95 % SDS-PAGE. ab92769 is at least 98 percent pure as determined by SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Escherichia coli
製品の詳細	Recombinant <i>E. coli</i> Ecotin protein

特性

Our **Abpromise guarantee** covers the use of **ab92769** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	Kinetic Data: 1:1 inhibition of trypsin (mole per mole) Previously labelled as E. coli serine protease inhibitor.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 3 Constituent: 0.0036% Hydrochloric acid This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

関連情報

関連性	E. coli serine protease inhibitor (Ecotin) inhibits pancreatic serine proteases: inhibits chymotrypsin, trypsin, elastases, factor X, kallikrein as well as a variety of other proteases. The
-----	---

strength of inhibition does not appear to be correlated with a particular protease specificity.

細胞内局在

Periplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors