

Recombinant Hepatitis B virus Hepatitis B Surface Antigen AD protein ab193473

1 References

製品の詳細

製品名	Recombinant Hepatitis B virus Hepatitis B Surface Antigen AD protein
精製度	> 95 % SDS-PAGE. Recombinant Hepatitis B surface antigen AD expressed in <i>Saccharomyces cerevisiae</i> and purified via ionic exchange chromatography, size exclusion chromatography and sterile filtration.
発現系	<i>Saccharomyces cerevisiae</i>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Hepatitis B virus
予測される分子量	24 kDa
配列の追加情報	Surface antigen Ad.

特性

Our **Abpromise guarantee** covers the use of **ab193473** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE ELISA Western blot
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C. Do Not Freeze. Preservative: 0.09% Sodium azide Constituent: 99% PBS
------------	---

関連情報

関連性	Hepatitis B Virus (HBV) infection induces a disease state which manifests itself in a variety of ways, characterized by the extent of liver damage, inflammation and viral persistence. HBV infection is also associated with a 100 fold increased risk of hepatocellular carcinoma and currently infects over 250 million people worldwide. HBV has a partially double stranded 3.2 kilobase DNA genome which contains four open reading frames. One of these encodes a 154 amino acid protein called the HBx protein. HBx has been shown to be a transcriptional transactivator of both viral and cellular promoters. Lacking a DNA binding domain and nuclear localization signal, HBx is believed to exert transcriptional activity through protein protein interaction.
細胞内局在	Cell Membrane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors